

Biotronix
Healthcare®

RESPIRATORY CATALOG



Biotronix Healthcare Industries Inc.
Ft. Lauderdale, Florida USA



CONTENT

Aerosol Therapy 4

Aerosol Masks

Aerosol Nebulizer Kits

Oxygen Therapy 10

Airways

Oxygen Masks

Nasal Cannulas

Oxygen Supply Tubing

Oxygen Corrugated Tubing

Advanced Breathing Circuit

Coaxial Breathing Circuits

Circuits 29

Advanced Breathing Circuit

Coaxial Breathing Circuits

Expandable Breathing Circuits

Corrugated Breathing Circuits

Biotronix **Healthcare®**

Aerosol Masks

Biotronix Healthcare offers a full range of aerosol masks that maximizes patient and clinical outcomes. They are designed for use with Biotronix Healthcare nebulizers and are available in adult, pediatric and infant sizes.


Standard Aerosol Mask

- Designed for patient comfort
- Clear, soft vinyl
- Extremely lightweight and flexible
- Adjustable nose clip
- Accepts 22mm ID aerosol tubing



Biotronix Healthcare Standard Aerosol Mask

REF #	DESCRIPTION	QUANTITY (per box)	Quantity (per case)
RE334001	Adult Standard Aerosol Mask Only	N/A	50/Case
RE334003	Pediatric Standard Aerosol Mask Only	N/A	50/Case
RE334000	Infant Standard Aerosol Mask Only	N/A	50/Case



Aerosol Mask-Elongated

- Under the chin designed for patient comfort
- Clear, soft vinyl
- Extremely lightweight and flexible
- Adjustable nose clip
- Accepts 22mm ID aerosol tubing



Biotronix Healthcare Aerosol Mask-Elongated

REF #	DESCRIPTION	QUANTITY (per box)	Quantity (per case)
RE334002	Adult Elongated Aerosol Mask only	N/A	50/Case
RE334004	Pediatric Elongated Aerosol Mask only	N/A	50/Case
RE334005	Infant Elongated Aerosol Mask only	N/A	50/Case

Face Tent

- Clear, soft vinyl
- Elastic strap
- Under the chin design
- For high humidity aerosol therapy
- Thick wall for mask stability



Biotronix Healthcare Face Tent

REF #	DESCRIPTION	QUANTITY (per box)	Quantity (per case)
RE336001	Face Tent, Under the Chin, no tubing included	N/A	50/Case
RE336002	Face Tent, Under the Chin, with 6 Ft Corrugated Tubing	N/A	50/Case
RE335003	Face Tent, Under the Chin, with 6 inch resevoir tube, 7 FT No Crush OxygenSupply Tubing and O2 adapter	N/A	50/Case


Tracheostomy Mask

- Clear, soft vinyl
- Designed for patient comfort
- Helix connector for 360° swivel
- 22mm outer connector



Biotronix Healthcare Tracheostomy Mask

REF #	DESCRIPTION	QUANTITY (per box)	Quantity (per case)
RE335001	Adult Tracheostomy Mask Only	N/A	50/Case
RE335002	Pediatric Tracheostomy Mask Only	N/A	50/Case
RE33503	Infant Tracheostomy Mask Only	N/A	50/Case



Biotronix
Healthcare®
Aerosol Masks

Biotronix Healthcare®

Aerosol Nebulizer Kits

Biotronix Healthcare offers a complete range of Nebulizer kits that provides excellent patient outcome and assists healthcare facilities efficiency. Our products include nebulizers with multiple hose, connector and volume configurations as well as configurations which help patients experience a comfortable therapy session.



Nebulizer with Mouth Piece

- Designed for patient comfort
- Clear, soft vinyl
- Flow Rate 4-8/ lpm
- Extremely lightweight and flexible
- Adjustable nose clip
- Accepts 22mm ID aerosol tubing
- Easy seal 6ml capacity spill proof jar
- Available with Standard and Universal Connector



Biotronix Healthcare Nebulizers with Mouth Piece

REF #	DESCRIPTION	Quantity (per case)
RE334001S	Nebulize with T Piece, Mouth Piece, & 7ft Tubing with Standard connector	50/Case
RE334001U	Nebulize with T Piece, Mouth Piece, & 7ft Tubing with Universal connector	50/Case
RE334002-S	Nebulize with T Piece, Mouth Piece, Reservoir Tube, & 7ft Tubing with Standard connector	50/Case
RE334003	Nebulizer with T Piece, & Mouth Piece	50/Case

Nebulizer Kit with Aerosol Mask

- Under the chin designed for patient comfort
- Clear, soft vinyl
- Flow Rate 4-8/ lpm
- Extremely lightweight and flexible
- Adjustable nose clip
- Accepts 22mm ID aerosol tubing
- Easy Seal 6ml capacity spill proof jar
- Available with Standard and Universal Connector




Biotronix Healthcare Nebulizers with Aerosol Masks

REF #	DESCRIPTION	Quantity (per case)
RE343001S	Nebulizer with Adult Standard Mask & 7ft No Crush Tubing with OxyTru Standard Connector	50/Case
RE343001SS	Nebulizer with Adult Standard Mask w/Swivel Connector & 7ft No Crush Tubing with OxyTru Standard Connector	50/Case
RE343002-S	Nebulizer with Adult Elongated Mask & 7ft No Crush Tubing with OxyTru Standard Connector	50/Case
RE343002-SS	Nebulizer with Adult Elongated Mask w/Swivel Connector & 7ft No Crush Tubing with OxyTru Standard Connector	50/Case
RE343002-U	Nebulizer with Adult Elongated Mask & 7ft No Crush Tubing with OxyOne Universal Connector	50/Case
RE343002-US	Nebulizer with Adult Elongated Mask w/Swivel Connector & 7ft No Crush Tubing with OxyOne Universal Connector	50/Case
RE343003-S	Nebulizer with Pediatric Standard Mask & 7ft No Crush Tubing with OxyTru Standard Connector	50/Case
RE343003-SS	Nebulizer with Pediatric Standard Mask w/Swivel Connector & 7ft No Crush Tubing with OxyTru Standard Connector	50/Case
RE343004-S	Nebulizer with Pediatric Elongated Mask & 7ft No Crush Tubing with OxyTru Standard Connector	50/Case
RE343004-SS	Nebulizer with Pediatric Elongated Mask w/Swivel Connector & 7ft No Crush Tubing with OxyTru Standard Connector	50/Case
RE343004-U	Nebulizer with Pediatric Elongated Mask & 7ft No Crush Tubing with OxyOne Universal Connector	50/Case
RE343004-US	Nebulizer with Pediatric Elongated Mask w/Swivel Connector & 7ft No Crush Tubing with OxyOne Universal Connector	50/Case

Biotronix Healthcare®

Aerosol Nebulizer Kits



Biotronix Healthcare®

Airways

Nasopharyngeal Airways

- Made from soft PVC
- Rounded tip for gentle insertion
- Unique design for secure placement
- Flexible for excellent patient comfort
- Available in a full range of sizes including pediatric.
- Latex Free



Biotronix Healthcare Nasopharyngeal Airways Trumpet Style

REF #	SIZE/ID (FR)	SIZE/ID (MM)	QUANTITY (per box)	Quantity (per case)
AN321030	12	3.0mm	10/Box	10 Boxes/Case
AN321035	14	3.5mm	10/Box	10 Boxes/Case
AN321040	16	4.0mm	10/Box	10 Boxes/Case
AN321045	18	4.5mm	10/Box	10 Boxes/Case
AN321050	20	5.0mm	10/Box	10 Boxes/Case
AN321055	22	5.5mm	10/Box	10 Boxes/Case
AN321060	24	6.0mm	10/Box	10 Boxes/Case
AN321065	26	6.5mm	10/Box	10 Boxes/Case
AN321070	28	7.0mm	10/Box	10 Boxes/Case
AN321075	30	7.5mm	10/Box	10 Boxes/Case
AN321080	32	8.0mm	10/Box	10 Boxes/Case
AN321085	34	8.5mm	10/Box	10 Boxes/Case
AN321090	36	9.0mm	10/Box	10 Boxes/Case

Nasopharyngeal Airway w/Adjustable Flange

- Semi-rigid polyethylene construction
- Two channel
- Individually packaged
- Non-Toxic
- Latex Free



Biotronix Healthcare Nasopharyngeal Airway w/Adjustable Flange

REF #	SIZE/ID (FR)	SIZE/ID (MM)	QUANTITY. (per box)	Quantity (per case)
AN322030	12	3.0mm	10/Box	10 Boxes/Case
AN322035	14	3.5mm	10/Box	10 Boxes/Case
AN322040	16	4.0mm	10/Box	10 Boxes/Case
AN322045	18	4.5mm	10/Box	10 Boxes/Case
AN322050	20	5.0mm	10/Box	10 Boxes/Case
AN322055	22	5.5mm	10/Box	10 Boxes/Case
AN322060	24	6.0mm	10/Box	10 Boxes/Case
AN322065	26	6.5mm	10/Box	10 Boxes/Case
AN322070	28	7.0mm	10/Box	10 Boxes/Case
AN322075	30	7.5mm	10/Box	10 Boxes/Case
AN322080	32	8.0mm	10/Box	10 Boxes/Case
AN322085	34	8.5mm	10/Box	10 Boxes/Case
AN322090	36	9.0mm	10/Box	10 Boxes/Case

Guedell Airways

- Nontoxic polyethylene construction
- Semi rigid, enclosed centre channel
- Color coded
- Individually packaged
- Latex Free



Biotronix Healthcare Guedel Airway

REF #	SIZE	LENGTH	COLOR	QUANTITY (per box)	Quantity (per case)
AN3110	00	40mm	PINK	10/Box	10 Boxes/Case
AN3110	0	50mm	BLUE	10/Box	10 Boxes/Case
AN3110	1	60mm	BLACK	10/Box	10 Boxes/Case
AN3110	2	70mm	WHITE	10/Box	10 Boxes/Case
AN3110	3	80mm	GREEN	10/Box	10 Boxes/Case
AN3110	4	90mm	YELLOW	10/Box	10 Boxes/Case
AN3110	5	100mm	RED	10/Box	10 Boxes/Case
AN3110	6	110mm	ORANGE	10/Box	10 Boxes/Case

Berman Airways

- Semi-rigid polyethylene construction
- Two channel
- Individually packaged
- Non-Toxic
- Latex Free



Biotronix Healthcare Berman Airway

REF #	SIZE	LENGTH	QUANTITY (per box)	Quantity (per case)
AN3120	00	40mm	10/Box	10 Boxes/Case
AN3120	0	50mm	10/Box	10 Boxes/Case
AN3120	1	60mm	10/Box	10 Boxes/Case
AN3120	2	70mm	10/Box	10 Boxes/Case
AN3120	3	80mm	10/Box	10 Boxes/Case
AN3120	4	90mm	10/Box	10 Boxes/Case
AN3120	5	100mm	10/Box	10 Boxes/Case
AN3120	6	110mm	10/Box	10 Boxes/Case

Biotronix Healthcare®

Oxygen Masks

Biotronix Healthcare has designed a complete range of oxygen masks to fulfill any clinical need. Our masks are soft, clear, light, odorless, made from high quality resins; the adjustable mask straps and nose clip create a better and more comfortable fit for the patient. Together with our innovative no-crush/non-kink oxygen tubing we created oxygen masks to meet all your clinical needs.



OxyTru Standard

OxyOne Universal

High Concentration Oxygen Masks

- Clear, soft vinyl for patient comfort and visual assessment
- Adaptor swivels to accommodate patient position
- Adjustable nose clip assures comfortable fit
- 750-mL reservoir bag
- Complete with 7-ft standard oxygen supply tubing or No crush oxygen tubing
- With Standard or Universal connectors



Biotronix Healthcare High Concentration Oxygen Masks

REF #	DESCRIPTION	Quantity (per case)
RES361007S	Adult High-Concentration Standard Oxygen Mask, with Reservoir Bag, 7ft No Crush Oxygen Supply Tubing, Standard Connector	50 Boxes/Case
RES361008S	Adult High-Concentration Elongated Oxygen Mask, with Reservoir Bag, 7ft No Crush Oxygen Supply Tubing, Standard Connector	50 Boxes/Case
RES361007U	Adult High-Concentration Standard Oxygen Mask, with Reservoir Bag, 7ft No Crush Oxygen Supply Tubing, Universal Connector	50 Boxes/Case
RES361008U	Adult High-Concentration Elongated Oxygen Mask, with Reservoir Bag, 7ft No Crush Oxygen Supply Tubing, Universal Connector	50 Boxes/Case
RES361009S	Pediatric High-Concentration Standard Oxygen Mask, with Reservoir Bag, 7ft No Crush Oxygen Supply Tubing, Standard Connector	50 Boxes/Case
RES361010S	Pediatric High-Concentration Elongated Oxygen Mask, with Reservoir Bag, 7ft No Crush Oxygen Supply Tubing, Standard Connector	50 Boxes/Case
RES361009U	Pediatric High-Concentration Standard Oxygen Mask, with Reservoir Bag, 7ft No Crush Oxygen Supply Tubing, Universal Connector	50 Boxes/Case
RES361010U	Pediatric High-Concentration Elongated Oxygen Mask, with Reservoir Bag, 7ft No Crush Oxygen Supply Tubing, Universal Connector	50 Boxes/Case
RES361011S	Infant High-Concentration Standard Oxygen Mask, with Reservoir Bag, 7ft No Crush Oxygen Supply Tubing, Standard Connector	50 Boxes/Case
RES361012S	Infant High-Concentration Elongated Oxygen Mask, with Reservoir Bag, 7ft No Crush Oxygen Supply Tubing, Standard Connector	50 Boxes/Case
RES361011U	Infant High-Concentration Standard Oxygen Mask, with Reservoir Bag, 7ft No Crush Oxygen Supply Tubing, Universal Connector	50 Boxes/Case
RES361012U	Infant High-Concentration Elongated Oxygen Mask, with Reservoir Bag, 7ft No Crush Oxygen Supply Tubing, Universal Connector	50 Boxes/Case

Medium Concentration Oxygen Masks

- Standard and under-the-chin designs for excellent fit on wide range of face sizes
- Clear, soft vinyl for patient comfort and visual assessment
- Adjustable nose clip assures comfortable fit
- Complete with 7-ft standard oxygen supply tubing or No crush oxygen tubing
- With Standard or Universal connectors



Biotronix Healthcare Medium Concentration Oxygen Masks

REF #	DESCRIPTION	Quantity (per case)
RES331001S	Adult Medium-Concentration Standard Oxygen Mask, with 7ft No Crush Tubing, Standard Connector	50 Boxes/Case
RES331001U	Adult Medium-Concentration Standard Oxygen Mask, with 7ft No Crush Tubing, Universal Connector	50 Boxes/Case
RES331003S	Adult Medium-Concentration Elongated Oxygen Mask, with 7ft No Crush Tubing, Standard Connector	50 Boxes/Case
RES331003U	Adult Medium-Concentration Elongated Oxygen Mask, with 7ft No Crush Tubing, Universal Connector	50 Boxes/Case
RES331002	Adult Medium-Concentration Standard Oxygen Mask, No Tubing	50 Boxes/Case
RES331004	Adult Medium-Concentration Elongated Oxygen Mask, No Tubing	50 Boxes/Case
RES33105S	Pediatric Medium-Concentration Standard Oxygen Mask, with 7ft No Crush Tubing, Standard Connector	50 Boxes/Case
RES331005U	Pediatric Medium-Concentration Standard Oxygen Mask, with 7ft No Crush Tubing, Universal Connector	50 Boxes/Case
RES331007S	Pediatric Medium-Concentration Elongated Oxygen Mask, with 7ft No Crush Tubing, Standard Connector	50 Boxes/Case
RES331007U	Pediatric Medium-Concentration Elongated Oxygen Mask, with 7ft No Crush Tubing, Universal Connector	50 Boxes/Case
RES331008	Pediatric Medium-Concentration Standard Oxygen Mask, No Tubing	50 Boxes/Case
RES331009	Pediatric Medium-Concentration Elongated Oxygen Mask, No Tubing	50 Boxes/Case
RES331010	Infant Medium-Concentration Oxygen Mask, No Tubing	50 Boxes/Case
RES331005IS	Infant Medium-Concentration Oxygen Mask, with 7ft No Crush Tubing, Standard Connector	50 Boxes/Case

Non-Rebreathing Oxygen Masks

- Clear, soft vinyl for patient comfort and visual assessment
- Low-resistance check valve does not contain natural rubber latex prevents rebreathing and allows exhaled gas to escape
- Adjustable nose clip assures comfortable fit
- 750-mL reservoir bag
- Complete with 7-ft standard oxygen supply tubing or No crush oxygen tubing



Biotronix Healthcare Non-Rebreathing Oxygen Masks

REF #	DESCRIPTION	Quantity (per case)
RES361007S	Adult Non-Rebreathing Oxygen Mask, 7ft No Crush Tubing, Standard Connector	50 Boxes/Case
RES361007U	Adult Non-Rebreathing Oxygen Mask, 7ft No Crush Tubing, Universal Connector	50 Boxes/Case
RES361014S	Pediatric Non-Rebreathing Oxygen Mask, 7ft No Crush Tubing, Standard Connector	50 Boxes/Case
RES361014U	Pediatric Non-Rebreathing Oxygen Mask, 7ft No Crush Tubing, Universal Connector	50 Boxes/Case
RES361025S	Infant Non-Rebreathing Oxygen Mask, 7ft No Crush Tubing, Standard Connector	50 Boxes/Case
RES361025U	Infant Non-Rebreathing Oxygen Mask, 7ft No Crush Tubing, Universal Connector	50 Boxes/Case

Venturi Oxygen Masks

- Clear, soft vinyl for patient comfort and visual assessment
- Features 6 color-coded diluters: 24%, 28%, 31%, 35%, 40% and 50%
- Adjustable nose clip assures comfortable fit
- Complete with 7-ft standard oxygen supply tubing or No crush oxygen tubing
- With Standard or Universal connectors



Biotronix Healthcare Venturi Oxygen Masks

REF #	DESCRIPTION	Quantity (per case)
RES332001S	Adult Venturi Oxygen Mask, with 6 Color Coded Diluters 24%, 28%, 31%, 35%, 40% and 50%, adaptor for high-humidity, 7ft No Crush Oxygen Tubing, Standard Connector	50 Boxes/Case
RES332001U	Adult Venturi Oxygen Mask, with 6 Color Coded Diluters 24%, 28%, 31%, 35%, 40% and 50%, adaptor for high-humidity, 7ft No Crush Oxygen Tubing, Universal Connector	50 Boxes/Case
RES332003U	Pediatric Venturi Oxygen Mask, with 6 Color Coded Diluters 24%, 28%, 31%, 35%, 40% and 50%, adaptor for high-humidity, 7ft No Crush Oxygen Tubing, Universal Connector	50 Boxes/Case

Multi-Vent Oxygen Masks

- Locking ring secures flow setting. Includes 2 color coded diluters: green for low-concentration, white for medium-concentration
- Clear, soft vinyl for patient comfort and visual assessment
- Complete with 7-ft standard oxygen supply tubing or No crush oxygen tubing
- With Standard or Universal connectors



Biotronix Healthcare Multi-Vent Oxygen Masks

REF #	DESCRIPTION	Quantity (per case)
RES332002S	Adult Multi-Vent Mask, Includes 2 color coded diluters: green for low-concentration, white for medium-concentration, 7ft No Crush Tubing, Standard Connector	50 Boxes/Case
RES332002U	Adult Multi-Vent Mask, Includes 2 color coded diluters: green for low-concentration, white for medium-concentration, 7ft No Crush Tubing, Universal Connector	50 Boxes/Case
RES332004S	Pediatric Multi-Vent Mask, Includes 2 color coded diluters: green for low-concentration, white for medium-concentration, 7ft No Crush Tubing, Universal Connector	50 Boxes/Case

Tracheostomy Masks

- Clear, soft vinyl for patient comfort and visual assessment
- Helix Swivel 22mm Connector
- Secure Strap



Biotraonix Healthcare Tracheostomy Masks

REF #	DESCRIPTION	Quantity (per case)
AN415001	Tracheostomy Mask, Adult,	50/Case
AN415006	Tracheostomy Mask, Pediatric,	50/Case

Biotronix **Healthcare®**

Nasal Cannulas

Biotronix Healthcare has an extensive range of Nasal Cannulas for all types of healthcare facilities as well as homecare use. Our Nasal Cannulas are manufactured using Non-Crush tubing with soft curved prongs for enhanced patient safety and comfort as well as standard tubing with standard prongs for a cost-effective solution.



OxyTru Standard

OxyOne Universal

SoftFlow II DEHP Free Nasal Cannula

- DEHP Free
- The SoftFlow extra soft nasal prongs designed for patient comfort
- Extremely lightweight and flexible
- Over the ear design





Biotronix Healthcare Soft-Flow II DEHP Free Nasal Cannula

REF #	DESCRIPTION	Quantity (per case)
RE356001S	Adult SoftFlow II Over-the-Ear Nasal Cannula, 7ft No Crush Tubing, Standard connector,	50/Case
RE356001U	Adult SoftFlow II Over-the-Ear Nasal Cannula, 7ft No Crush Tubing, Universal connector	50/Case
RE356001.14S	Adult SoftFlow II Over-the-Ear Nasal Cannula, 14ft No Crush Tubing, Standard connector,	50/Case
RE356002S	Pediatric SoftFlow II Over-the-Ear Nasal Cannula, 7ft No Crush Tubing, Standard connector,	50/Case
RE356002U	Pediatric SoftFlow II Over-the-Ear Nasal Cannula, 7ft No Crush Tubing, Universal connector	50/Case
RE356002.14S	Pediatric SoftFlow II Over-the-Ear Nasal Cannula, 14ft No Crush Tubing, Standard connector,	50/Case
RE356003S	Infant SoftFlow II Over-the-Ear Nasal Cannula, 7ft No Crush Tubing, Standard connector,	50/Case
RE356004s	Neonate SoftFlow II Over-the-Ear Nasal Cannula, 7ft No Crush Tubing, Standard connector,	50/Case

Soft Flow Nasal Cannula

- The SoftFlow extra soft nasal prongs designed for patient comfort
- Extremely lightweight and flexible
- Over the ear design
- Crush proof and kink resistant tubing



Biotronix Healthcare Soft Flow Nasal Cannula

REF #	DESCRIPTION	Quantity (per case)
RE354000S	Adult SoftFlow Over-the-Ear Nasal Cannula, 7ft Standard Tubing, Standard connector	50/Case
RE354001S	Adult SoftFlow Over-the-Ear Nasal Cannula, 7ft No Crush Tubing, Standard connector,	50/Case
RE354001U	Adult SoftFlow Over-the-Ear Nasal Cannula, 7ft No Crush Tubing, Universal connector	50/Case
RE354001.14S	Adult SoftFlow Over-the-Ear Nasal Cannula, 14ft No Crush Tubing, Standard connector,	50/Case
RE354001.14U	Adult SoftFlow Over-the-Ear Nasal Cannula, 14ft No Crush Tubing, Universal connector	50/Case
RE3564001.25S	Adult SoftFlow Over-the-Ear Nasal Cannula, 25ft No Crush Tubing, Standard connector,	25/Case
RE354001.25U	Adult SoftFlow Over-the-Ear Nasal Cannula, 25ft No Crush Tubing, Universal connector	25/Case
RE354001.50S	Adult SoftFlow Over-the-Ear Nasal Cannula, 50ft No Crush Tubing, Standard connector,	25/Case
RE354001.50U	Adult SoftFlow Over-the-Ear Nasal Cannula, 50ft No Crush Tubing, Universal connector	25/Case
RE354002S	Pediatric SoftFlow Over-the-Ear Nasal Cannula, 7ft No Crush Tubing, Standard connector,	50/Case
RE354002U	Pediatric SoftFlow Over-the-Ear Nasal Cannula, 7ft No Crush Tubing, Universal connector	50/Case
RE354002.14S	Pediatric SoftFlow Over-the-Ear Nasal Cannula, 14ft No Crush Tubing, Standard connector,	50/Case
RE354003S	Infant SoftFlow Over-the-Ear Nasal Cannula, 7ft No Crush Tubing, Standard connector,	50/Case

Over the Ear Nasal Cannulas

- Designed for patient comfort
- Nontoxic medical grade PVC
- Extremely lightweight and flexible
- Over the ear design

Biotronix Healthcare Over the Ear Nasal Cannulas

REF #	DESCRIPTION	Quantity (per case)
RE351000SU	Adult Over-the-Ear Nasal Cannula, 7ft Standard Tubing, Universal connector, Non-Flared Nasal Tip	50/Case
RE351000S	Adult Over-the-Ear Nasal Cannula, 7ft Standard Tubing, Standard connector, Non-Flared Nasal Tip	50/Case
RE351001S	Adult Over-the-Ear Nasal Cannula, 7ft No Crush Tubing, Standard connector, Non-Flared Nasal Tip	50/Case
RE351001F	Adult Over-the-Ear Nasal Cannula, 7ft No Crush Tubing, Standard connector, Flared Nasal Tip	50/Case
RE351001U	Adult Over-the-Ear Nasal Cannula, 7ft No Crush Tubing, Universal connector, Non-Flared Nasal Tip	50/Case
RE351001.14S	Adult Over-the-Ear Nasal Cannula, 14 ft No Crush Tubing, Standard connector, Non-Flared Nasal Tip	50/Case
RE351001.25S	Adult Over-the-Ear Nasal Cannula, 25 ft No Crush Tubing, Standard connector, Non-Flared Nasal Tip	50/Case
RE351001.50S	Adult Over-the-Ear Nasal Cannula, 50 ft No Crush Tubing, Standard connector, Non-Flared Nasal Tip	25/Case
RE351002S	Pediatric Over-the-Ear Nasal Cannula, 7ft No Crush Tubing, Standard connector, Non-Flared Nasal Tip	50/Case
RE351002U	Pediatric Over-the-Ear Nasal Cannula, 7ft No Crush Tubing, Universal connector, Non-Flared Nasal Tip	50/Case
RE351003S	Infant Over-the-Ear Nasal Cannula, 7ft No Crush Tubing, Standard connector, Non-Flared Nasal Tip	50/Case
RE351003U	Infant Over-the-Ear Nasal Cannula, 7ft No Crush Tubing, Universal connector, Non-Flared Nasal Tip	50/Case

Biotronix Healthcare®

Oxygen Supply Tubing

Biotronix Healthcare has an extensive range of oxygen tubing for all types of healthcare facilities as well as homecare use. No Crush tubing for enhanced patient safety and comfort as well as quality standard tubing. A variety of lengths and connectors types to meet patient needs and improving your healthcare facilities efficiency.




- Extremely Lightweight and Flexible
- No Crush Tubing for enhanced patient safety and Standard Smooth Bore Tubing
- Standard Cone Connector and Universal Connector that attached directly to flow meter connector
- Clear tubing or Green tubing for patient safety

Biotronix Healthcare Oxygen Supply Tubing

REF #	STANDARD TUBING CLEAR	Quantity (per case)
RE361007S	7ft No-Crush Clear Tubing, Standard Connector	50/Case
RE362007S	7ft Standard Clear Tubing, Standard Connector	50/Case
RE362007U	7ft Standard Clear Tubing, Universal Connector	50/Case
RE361014S	14ft No-Crush Clear Tubing, Standard Connector	50/Case
RE363014S	14ft No-Crush Green Tubing, Standard Connector	50/Case
RE361025S	25ft No-Crush Clear Tubing, Standard Connector	25/Case
RE361025U	25ft No Crush Clear Tubing, Universal Connector	25/Case
RE363025S	25ft No-Crush Green Tubing, Standard Connector	25/Case
RE362025S	25ft Standard Clear Tubing, Standard Connector	25/Case
RE361030S	30ft No Crush Clear Tubing, Standard Connector	20/Case
RE361040S	40ft No-Crush Clear Tubing, Standard Connector	15/Case
RE363040S	40ft No-Crush Green Tubing, Standard Connector	15/Case
RE361050S	50ft No-Crush Clear Tubing, Standard Connector	25/Case
RE363050S	50ft No-Crush Green Tubing, Standard Connector	25/Case
RE362050S	50ft Standard Clear Tubing, Standard Connector	25/Case

Biotronix Healthcare®

Oxygen Corrugated Tubing



The Biotronix Healthcare Oxygen Corrugated Tubing is disposable, flexible, light, available in precut lengths and segmented every 6 inches for easy cutting. Our corrugated tubing is perfect for aerosol / IPPB applications offering flexibility and convenience to healthcare facilities.

- Tapered Cuffs
- Segmented every 6"
- Dispenser Box
- 22mm attachment
- Individually packaged

Biotronix Healthcare Oxygen Corrugated Tubing

REF #	STANDARD TUBING CLEAR	Quantity (per case)
RE37001	Corrugated Oxygen Supply Tubing 100 Feet, Cuttable Cuff, Segmented every 6"	1/Case
RE37002	Corrugated Oxygen Supply Tubing 72", Cuttable Cuff, Segmented every 6"	50/Case
RE37003	Corrugated Oxygen Supply Tubing 60", Cuttable Cuff, Segmented every 6"	50/Case
RE37004	Corrugated Oxygen Supply Tubing 48", Cuttable Cuff, Segmented every 6"	50/Case
RE37005	Corrugated Oxygen Supply Tubing 6", No Cuttable	50/Case



Biotronix[®]
Healthcare

BI-LUMEN DESIGN
ADVANCED
BREATHING CIRCUIT



Clinical Benefits:

1) The VeroX circuit is ideal for both Adults and Pediatrics:

Since the total circuit volume is much lower than a standard Adult corrugated circuit (either circle or coaxial), the VeroX circuits are ideally suited for use with both Adults and Pediatrics, which then avoids the need to stock two different circuits and saves time by eliminating the need to change circuits between these two patient groups.

2) Ideal for Low/Ultra Low flow Anesthesia:

a) Quicker wash-in/wash-out times: - with a much lower total circuit volume than a standard corrugated circle circuit, this allows for quicker wash in/wash out times.

b) Reduces rain out - One of the side effects of running lower flows is the potential for rain-out in the circuit. Cold inspiratory gases are warmed through the inner wall by the returned gases, which helps reduce the overall circuit temperature differential to room air temperature, which in turn reduces rain-out.

3) The VeroX is lighter, smaller and creates a less cluttered work area:

Approximately half the weight of many corrugated Adult circle circuits and is even smaller and lighter than a Coaxial circuit. This means that it is easier to place on the patient, reduces torque on the ET Tube/LMA and also provides a cleaner work environment for the Anesthesiologist.

4) Guaranteed Quality through 100% leak testing:

Unlike most manufacturers, every circuit is fully leak tested as part of our stringent quality control measures. Our leak test acceptance levels are far in excess of the current ISO 5367

standards, which gives the user complete peace of mind that the circuit will pass the most advanced anesthetic machine tests for compliance and leaks and will provide high performance under the most extreme clinical situations such as ultra low flow anesthesia and pediatric anesthesia.

5) VeroX circuit uses Silver Ion technology which helps in the reduction of many different types of bacteria found in hospitals in as little as 15 minutes and be 99.99% effective within 2 hours.

6) The Integrated Gas sampling lines positioned in the centre of the tidal flow, which is where the best flow is found when tidal volumes are minimal. This will provide a better response for low respiratory output patients, compared to using the sample port on a filter.

Biotronix Healthcare VeroX Adult Bi-Lumen Breathing Circuits

REF #	TUBING	WYE	LENGTH	FITTINGS	FILTER/HMEF	LATEX FREE BREATHING BAG	GAS SAMPLING LINE	MASK	CASE QTY
ANOX71001	VeroX Bi-Lumen	Straight	60" Inch	Gas Sampling Elbow	-	-	-	-	20/Case
ANOX71002	VeroX Bi-Lumen	Straight	60" Inch	Gas Sampling Elbow	Bacterial Viral (BV)	-	-	-	20/Case
ANOX71003	VeroX Bi-Lumen	Straight	60" Inch	Gas Sampling Elbow	HMEF	-	-	-	20/Case
ANOX71004	VeroX Bi-Lumen	Straight	60" Inch	Gas Sampling Elbow	HMEF	3L	-	-	20/Case
ANOX71005	VeroX Bi-Lumen	Straight	60" Inch	Gas Sampling Elbow	HMEF	3L	10' M/M	Adult #5	20/Case
ANOX71005 P	VeroX Bi-Lumen	Straight	60" Inch	Gas Sampling Elbow	HMEF	3L	10' M/M	Ped #3	20/Case
ANOX71006	VeroX Bi-Lumen	Straight	60" Inch	Gas Sampling Elbow	HMEF	3L	10' M/M	-	20/Case
ANOX71007	VeroX Bi-Lumen	Straight	60" Inch	Gas Sampling Elbow	HMEF	3L	-	Adult #5	20/Case
ANOX71007 P	VeroX Bi-Lumen	Straight	60" Inch	Gas Sampling Elbow	HMEF	3L	-	Ped #3	20/Case
ANOX71008	VeroX Bi-Lumen	Straight	72" Inch	Gas Sampling Elbow	-	-	-	-	20/Case
ANOX71009	VeroX Bi-Lumen	Straight	72" Inch	Gas Sampling Elbow	Bacterial Viral (BV)	-	-	-	20/Case
ANOX71010	VeroX Bi-Lumen	Straight	72" Inch	Gas Sampling Elbow	HMEF	-	-	-	20/Case
ANOX71011	VeroX Bi-Lumen	Straight	72" Inch	Gas Sampling Elbow	HMEF	3L	-	-	20/Case
ANOX71012	VeroX Bi-Lumen	Straight	72" Inch	Gas Sampling Elbow	HMEF	3L	10' M/M	-	20/Case
ANOX71013	VeroX Bi-Lumen	Straight	72" Inch	Gas Sampling Elbow	Bacterial Viral (BV)	2L	10' M/M	Adult #6	20/Case
ANOX71014	VeroX Bi-Lumen	Straight	72" Inch	Gas Sampling Elbow	HMEF	3L	-	Adult #5	20/Case
ANOX71015	VeroX Bi-Lumen	Straight	72" Inch	Gas Sampling Elbow	Bacterial Viral (BV)	3L	10' M/M	Adult #5	20/Case
ANOX71016	VeroX Bi-Lumen	Straight	72" Inch	Gas Sampling Elbow	HMEF	3L	10' M/M	Adult #5	20/Case
ANOX71017	VeroX Bi-Lumen	Straight	72" Inch	Gas Sampling Elbow	HMEF	3L	10' M/M	Adult #6	20/Case
ANOX71018	VeroX Bi-Lumen	Straight	90" Inch	Gas Sampling Elbow	-	-	-	-	20/Case
ANOX71019	VeroX Bi-Lumen	Straight	90" Inch	Gas Sampling Elbow	Bacterial Viral (BV)	-	-	-	20/Case
ANOX71020	VeroX Bi-Lumen	Straight	90" Inch	Gas Sampling Elbow	HMEF	-	-	-	20/Case
ANOX71021	VeroX Bi-Lumen	Straight	90" Inch	Gas Sampling Elbow	HMEF	3L	-	-	20/Case
ANOX71022	VeroX Bi-Lumen	Straight	90" Inch	Gas Sampling Elbow	HMEF	3L	-	Adult #5	20/Case
ANOX71023	VeroX Bi-Lumen	Straight	90" Inch	Gas Sampling Elbow	HMEF	3L	10' M/M	Adult #5	20/Case



Biotronix
Healthcare®

Verox Bi-Lumen Circuit

Biotronix Healthcare®

Coaxial Breathing Circuits

Biotronix Healthcare's Coaxial Breathing Circuits are designed to prevent common problems like added dead space and weight of a circle circuit system. The tube within a tube design eliminates tangling between the expiratory and inspiratory limbs. The design reduces clutter and keeps order in the ICU and anesthesia field and provides heat exchange between the limbs. Multiple configurations are available with choices that including masks, BV filters/HMEF, breathing bags, gas sampling lines, and different lengths.



Adult Circuit Features & Benefits

- Lightweight and low dead space,
- The inspiratory and expiratory limb tube and tube design eliminates twisting.
- Produces heat and moisture conservation properties better than standard dual limb design.
- Built in gas sampling port.
- 100% leak tested
- Biotronix Healthcare Breathing bags are made from a synthetic latex free material
- Multiple configurations available including masks, gas sampling lines, length of circuit, filters, HMEF's, breathing bags and extensions.


Biotronix Healthcare Adult Coaxile Breathing Circuits

REF #	TUBING	WYE	LENGTH	FITTINGS	FILTER/HMEF	LATEX FREE BREATHING BAG	GAS SAMPLING LINE	MASK	CASE QTY
ANK51001	Coaxile	Straight	60" Inch	Gas Sampling Elbow	-	-	-	-	20/Case
ANK51002	Coaxile	Straight	60" Inch	Gas Sampling Elbow	Bacterial Viral (BV)	-	-	-	20/Case
ANK51003	Coaxile	Straight	60" Inch	Gas Sampling Elbow	HMEF	-	-	-	20/Case
ANK51004	Coaxile	Straight	60" Inch	Gas Sampling Elbow	HMEF	3L	-	-	20/Case
ANK51005	Coaxile	Straight	60" Inch	Gas Sampling Elbow	HMEF	3L	10' M/M	Adult #5	20/Case
ANK51006	Coaxile	Straight	60" Inch	Gas Sampling Elbow	HMEF	3L	10' M/M	-	20/Case
ANK51007	Coaxile	Straight	60" Inch	Gas Sampling Elbow	HMEF	3l	-	Adult #5	20/Case
ANK51008	Coaxile	Straight	72" Inch	Gas Sampling Elbow	-	-	-	-	20/Case
ANK51009	Coaxile	Straight	72" Inch	Gas Sampling Elbow	Bacterial Viral (BV)	-	-	-	20/Case
ANK51010	Coaxile	Straight	72" Inch	Gas Sampling Elbow	HMEF	-	-	-	20/Case
ANK51011	Coaxile	Straight	72" Inch	Gas Sampling Elbow	HMEF	3L	-	-	20/Case
ANK51012	Coaxile	Straight	72" Inch	Gas Sampling Elbow	HMEF	3L	10' M/M	-	20/Case
ANK51013	Coaxile	Straight	72" Inch	Gas Sampling Elbow	Bacterial Viral (BV)	2L	10' M/M	Adult #6	20/Case
ANK51014	Coaxile	Straight	72" Inch	Gas Sampling Elbow	HMEF	3L	-	Adult #5	20/Case
ANK51015	Coaxile	Straight	72" Inch	Gas Sampling Elbow	Bacterial Viral (BV)	3L	10' M/M	Adult #5	20/Case
ANK51016	Coaxile	Straight	72" Inch	Gas Sampling Elbow	HMEF	3L	10' M/M	Adult #5	20/Case
ANK51017	Coaxile	Straight	72" Inch	Gas Sampling Elbow	HMEF	3L	10' M/M	Adult #6	20/Case
ANK51018	Coaxile	Straight	90" Inch	Gas Sampling Elbow	-	-	-	-	20/Case
ANK51019	Coaxile	Straight	90" Inch	Gas Sampling Elbow	Bacterial Viral (BV)	-	-	-	20/Case
ANK51020	Coaxile	Straight	90" Inch	Gas Sampling Elbow	HMEF	-	-	-	20/Case
ANK51021	Coaxile	Straight	90" Inch	Gas Sampling Elbow	HMEF	3L	-	-	20/Case
ANK51022	Coaxile	Straight	90" Inch	Gas Sampling Elbow	HMEF	3L	-	Adult #5	20/Case
ANK51023	Coaxile	Straight	90" Inch	Gas Sampling Elbow	HMEF	3L	10' M/M	Adult #5	20/Case

Pediatric Circuit

Features & Benefits

- Lightweight and low dead space,
- Smaller diameter tubing designed for pediatric use
- The inspiratory and expiratory limb tube and tube design eliminates twisting.
- Produces heat and moisture conservation properties better than standard dual limb design.
- Built in gas sampling port.
- 100% leak tested
- Biotronix Healthcare Breathing bags are made from a synthetic latex free material
- Multiple configurations available including masks, gas sampling lines, length of circuit, filters, HMEF's, breathing bags and extensions.






Biotronix Healthcare Pediatric Coaxile Breathing Circuits

REF #	TUBING	WYE	LENGTH	FITTINGS	FILTER/HMEF	LATEX FREE BREATHING BAG	GAS SAMPLING LINE	MASK	CASE QTY
ANK51050	Coaxile	Straight	40" Inch	Gas Sampling Elbow	-	-	-	-	20/Case
ANK51051	Coaxile	Straight	40" Inch	Gas Sampling Elbow	Bacterial Viral (BV)	-	-	-	20/Case
ANK51052	Coaxile	Straight	40" Inch	Gas Sampling Elbow	HMEF	-	-	-	20/Case
ANK51053	Coaxile	Straight	40" Inch	Gas Sampling Elbow	Bacterial Viral (BV)	1L	10' M/M	Ped #3	20/Case
ANK51054	Coaxile	Straight	40" Inch	Gas Sampling Elbow	HMEF	1L	10' M/M	Ped #3	20/Case
ANK51055	Coaxile	Straight	40" Inch	Gas Sampling Elbow	Bacterial Viral (BV)	1L	10' M/M	-	20/Case
ANK51056	Coaxile	Straight	40" Inch	Gas Sampling Elbow	Bacterial Viral (BV)	1L	-	Ped #3	20/Case
ANK51057	Coaxile	Straight	60" Inch	Gas Sampling Elbow	-	-	-	-	20/Case
ANK51058	Coaxile	Straight	60" Inch	Gas Sampling Elbow	Bacterial Viral (BV)	-	-	-	20/Case
ANK51059	Coaxile	Straight	60" Inch	Gas Sampling Elbow	HMEF	-	-	-	20/Case
ANK51060	Coaxile	Straight	60" Inch	Gas Sampling Elbow	Bacterial Viral (BV)	1L	10' M/M	Ped #3	20/Case
ANK51061	Coaxile	Straight	60" Inch	Gas Sampling Elbow	HMEF	1L	10' M/M	Ped #3	20/Case
ANK51062	Coaxile	Straight	60" Inch	Gas Sampling Elbow	Bacterial Viral (BV)	1L	10' M/M		20/Case
ANK51063	Coaxile	Straight	60" Inch	Gas Sampling Elbow	Bacterial Viral (BV)	1L	-	Ped #3	20/Case

Biotronix ***Healthcare®***

Expandable Breathing Circuits

Biotronix Healthcare's adult expandable breathing circuits are lightweight and flexible decreasing the possibility of kinking. Our 22mm expandable tubing allows gas to freely pass and allows the anesthesia machine and other equipment to be located away from the patient's head. Multiple configurations are available including choice of length of circuit, anesthesia masks, gas sampling lines, filters/HMEF's, water traps and breather bags to provide a complete cost-effective kit. Custom kits available.




Features & Benefits:

- 22 mm expandable tubing
- Lightweight and flexible
- Allows the anesthesia machine and other equipment to be located away from the patient's head

Adult Expandable Breathing Circuits

REF #	TUBING	WYE	LENGTH	FITTINGS	FILTER/HMEF	LATEX FREE BREATHING BAG	GAS SAMPLING LINE	MASK	CASE QTY
ANEX51001	Expandable	Parallel	72" Inch	Gas Sampling Elbow	-	-	-	-	20/Case
ANEX51002	Expandable	Parallel	72" Inch	Gas Sampling Elbow	Bacterial Viral (BV)	-	-	-	20/Case
ANEX51003	Expandable	Parallel	72" Inch	Gas Sampling Elbow	HMEF	-	-	-	20/Case
ANEX51004	Expandable	Parallel	72" Inch	Gas Sampling Elbow	Bacterial Viral (BV)	3L	-	-	20/Case
ANEX51005	Expandable	Parallel	72" Inch	Gas Sampling Elbow	Bacterial Viral (BV)	3L	10' M/M	Adult #5	20/Case
ANEX51006	Expandable	Parallel	72" Inch	Gas Sampling Elbow	HMEF	3L	10' M/M	Adult #5	20/Case
ANEX51007	Expandable	Parallel	72" Inch	Gas Sampling Elbow	Bacterial Viral (BV)	3L	-	Adult #5	20/Case
ANEX51008	Expandable	Parallel	72" Inch	Gas Sampling Elbow	Bacterial Viral (BV)	-	10" M/M	-	20/Case
ANEX51009	Expandable	Parallel	90" Inch	Gas Sampling Elbow	-	-	-	-	20/Case
ANEX51010	Expandable	Parallel	90" Inch	Gas Sampling Elbow	Bacterial Viral (BV)	-	-	-	20/Case
ANEX51011	Expandable	Parallel	90" Inch	Gas Sampling Elbow	HMEF	-	-	-	20/Case
ANEX51012	Expandable	Parallel	90" Inch	Gas Sampling Elbow	Bacterial Viral (BV)	3L	-	-	20/Case
ANEX51013	Expandable	Parallel	90" Inch	Gas Sampling Elbow	Bacterial Viral (BV)	3L	10' M/M	Adult #5	20/Case
ANEX51014	Expandable	Parallel	90" Inch	Gas Sampling Elbow	HMEF	3L	10' M/M	Adult #5	20/Case
ANEX51015	Expandable	Parallel	90" Inch	Gas Sampling Elbow	Bacterial Viral (BV)	3L	-	Adult #5	20/Case
ANEX51016	Expandable	Parallel	90" Inch	Gas Sampling Elbow	Bacterial Viral (BV)	-	10" M/M	-	20/Case
ANEX51017	Expandable	Parallel	108" Inch	Gas Sampling Elbow		3L			20/Case
ANEX51018	Expandable	Parallel	108" Inch	Gas Sampling Elbow	Bacterial Viral (BV)	3L			20/Case
ANEX51019	Expandable	Swivel	108" Inch	Gas Sampling Elbow	HMEF	3L			20/Case
ANEX51020	Expandable	Swivel	108" Inch	Gas Sampling Elbow					20/Case



Features & Benefits:

- Configurations include choice of mask, breather bag, gas sampling line, filters, HMEF's, and water traps
- Latex free breathing bags
- Able to manufacture breathing circuit lengths from 40" to 120"
- 20 each/case

Note: Multiple Configurations Available. Custom Kits Available.

Pediatric Expandable Breathing Circuits

REF #	TUBING	WYE	LENGTH	FITTINGS	FILTER/HMEF	LATEX FREE BREATHING BAG	GAS SAMPLING LINE	MASK	CASE QTY
ANK51048	Expandable	Parallel	60" Inch	Gas Sampling Elbow					
ANK51049	Expandable	Parallel	60" Inch	Gas Sampling Elbow	Bacterial Viral (BV)	1L			
ANK51050	Expandable	Parallel	72" Inch	Gas Sampling Elbow	-	-	-	-	20/Case
ANK51051	Expandable	Parallel	72" Inch	Gas Sampling Elbow	Bacterial Viral (BV)	-	-	-	20/Case
ANK51052	Expandable	Parallel	72" Inch	Gas Sampling Elbow	HMEF ANGLED	-	-	-	20/Case
ANK51053	Expandable	Parallel	72" Inch	Gas Sampling Elbow	Bacterial Viral (BV)	1L	10' M/M	Ped #3	20/Case
ANK51054	Expandable	Parallel	72" Inch	Gas Sampling Elbow	Bacterial Viral (BV)	1L	10' M/M		20/Case
ANK51055	Expandable	Parallel	90" Inch	Gas Sampling Elbow	-	-	-	-	20/Case
ANK51056	Expandable	Parallel	90" Inch	Gas Sampling Elbow	Bacterial Viral (BV)	-	-	-	20/Case
ANK51057	Expandable	Parallel	90" Inch	Gas Sampling Elbow	HMEF ANGLED	-	-	-	20/Case
ANK51058	Expandable	Parallel	90" Inch	Gas Sampling Elbow	Bacterial Viral (BV)	1L	10' M/M	Ped #3	20/Case
ANK51058	Expandable	Parallel	90" Inch	Gas Sampling Elbow	Bacterial Viral (BV)	1L	10' M/M		20/Case



Biotronix
Healthcare®

Expandable Breathing Circuits


Biotronix Healthcare®

Corrugated Breathing Circuits

Biotronix Healthcare's adult expandable breathing circuits are lightweight and flexible decreasing the possibility of kinking. Our 22mm expandable tubing allows gas to freely pass and allows the anesthesia machine and other equipment to be located away from the patient's head. Multiple configurations are available including choice of length of circuit, anesthesia masks, gas sampling lines, filters/HMEF's, water traps and breather bags to provide a complete cost-effective kit. Custom kits available.

Features & Benefits:

- 22 mm expandable tubing
- Lightweight and flexible
- Allows the anesthesia machine and other equipment to be located away from the patient's head
- Configurations include choice of mask, breather bag, gas sampling line, filters, HMEF's, and water traps
- Latex free breathing bags
- Available sizes of breathing circuit lengths from 40" to 120"
- 20 each/case



Biotronix Healthcare Adult Corrugated Breathing Circuits

REF #	TUBING	WYE	LENGT H	FITTINGS	Expiratory Lip Water Trap	Inspiratory Limb Water Trap	FILTER/ HMEF	LATEX FREE BREATHING BAG	GAS SAMPLIN G LINE	MASK	CASE QTY
ANCO51000	Corrugated	Paralle l	48" Inch	Gas Sampling Elbow	Y	Y	-	-	-	-	20/Case
ANCO51001	Corrugated	Paralle l	60" Inch	Gas Sampling Elbow	-	-	-	-	-	-	20/Case
ANCO51002	Corrugated	Paralle l	60" Inch	Gas Sampling Elbow	-	-	Bacterial Viral (BV)	-	-	-	20/Case
ANCO51003	Corrugated	Paralle l	60" Inch	Gas Sampling Elbow	-	-	HMEF	-	-	-	20/Case
ANCO51004	Corrugated	Paralle l	60" Inch	Gas Sampling Elbow	-	-	Bacterial Viral (BV)	3L	-	-	20/Case
ANCO51005	Corrugated	Paralle l	60" Inch	Gas Sampling Elbow	-	-	Bacterial Viral (BV)	3L	10' M/M	Adult #5	20/Case
ANCO51006	Corrugated	Paralle l	60" Inch	Gas Sampling Elbow	-	-	HMEF	3L	10' M/M	Adult #5	20/Case
ANCO51007	Corrugated	Paralle l	60" Inch	Gas Sampling Elbow	-	-	Bacterial Viral (BV)	3l	-	Adult #5	20/Case
ANCO51008	Corrugated	Paralle l	60" Inch	Gas Sampling Elbow	-	-	Bacterial Viral (BV)	-	10" M/M	-	20/Case
ANCO51009	Corrugated	Paralle l	60" Inch	Gas Sampling Elbow	Y	Y	-	-	-	-	20/Case
ANCO51010	Corrugated	Paralle l	72" Inch	Gas Sampling Elbow	Y	Y	-	-	10" M/M	-	20/Case
ANCO51013	Corrugated	Paralle l	90" Inch	Gas Sampling Elbow	-	-	-	-	-	-	20/Case
ANCO51014	Corrugated	Paralle l	90" Inch	Gas Sampling Elbow	-	-	Bacterial Viral (BV)	-	-	-	20/Case
ANCO51015	Corrugated	Paralle l	90" Inch	Gas Sampling Elbow	-	-	HMEF	-	-	-	20/Case
ANCO51016	Corrugated	Paralle l	90" Inch	Gas Sampling Elbow	-	-	Bacterial Viral (BV)	3L	-	-	20/Case
ANCO51017	Corrugated	Paralle l	90" Inch	Gas Sampling Elbow	-	-	Bacterial Viral (BV)	3L	10' M/M	Adult #5	20/Case
ANCO51018	Corrugated	Paralle l	90" Inch	Gas Sampling Elbow	-	-	HMEF	3L	10' M/M	Adult #5	20/Case
ANCO51019	Corrugated	Paralle l	90" Inch	Gas Sampling Elbow	-	-	Bacterial Viral (BV)	3L	-	Adult #5	20/Case
ANCO51020	Corrugated	Paralle l	90" Inch	Gas Sampling Elbow	-	-	Bacterial Viral (BV)	-	10" M/M	-	20/Case

***NOTE:** Custom kits are available. Please contact your sales representative.

Biotronix Healthcare Pediatric Corrugated Breathing Circuits


REF #	TUBING	WYE	LENGT H	FITTINGS	Expiratory Lip Water Trap	Inspiratory Limb Water Trap	FILTER/ HMEF	LATEX FREE BREATHIN G BAG	GAS SAMPLIN G LINE	MASK	CASE QTY
ANCP51048	Corrugated	Paralle l	60" Inch	Gas Sampling Elbow			-	-	-	-	20/Case
ANCP51049	Corrugated	Paralle l	60" Inch	Gas Sampling Elbow			Bacterial Viral (BV)	-	-	-	20/Case
ANCP51050	Corrugated	Paralle l	60" Inch	Gas Sampling Elbow			Bacterial Viral (BV)	1L	-	-	20/Case
ANCP51051	Corrugated	Paralle l	60" Inch	Gas Sampling Elbow			Bacterial Viral (BV)	1L	10' M/M	Ped #3	20/Case
ANCP51052	Corrugated	Paralle l	60" Inch	Gas Sampling Elbow			HMEF ANGLED	1L	10' M/M	Ped #3	20/Case
ANCP51053	Corrugated	Paralle l	60" Inch	Gas Sampling Elbow			Bacterial Viral (BV)	1L	10' M/M		20/Case
ANCP51053	Corrugated	Paralle l	60" Inch	Gas Sampling Elbow	Y	Y	-	-	-	-	20/Case
ANCP51054	Corrugated	Paralle l	90" Inch	Gas Sampling Elbow			-	-	-	-	20/Case
ANCP51055	Corrugated	Paralle l	90" Inch	Gas Sampling Elbow			Bacterial Viral (BV)	-	-	-	20/Case
ANCP51056	Corrugated	Paralle l	90" Inch	Gas Sampling Elbow			Bacterial Viral (BV)	1L	-	-	20/Case
ANCP51057	Corrugated	Paralle l	90" Inch	Gas Sampling Elbow			Bacterial Viral (BV)	1L	10' M/M	Ped #3	20/Case
ANCP51058	Corrugated	Paralle l	90" Inch	Gas Sampling Elbow			HMEF ANGLED	1L	10' M/M	Ped #3	20/Case
ANCP51059	Corrugated	Paralle l	90" Inch	Gas Sampling Elbow			Bacterial Viral (BV)	1L	10' M/M		20/Case

***NOTE:** Custom kits are available. Please contact your sales representative.



Biotronix
Healthcare®

Corrugated Breathing Circuits





Biotronik Healthcare

Biotronix Healthcare Industries Inc.
Ft. Lauderdale, Florida USA