


Biotronix
Healthcare®

ANESTHESIA

CATALOG



Biotronix Healthcare Industries Inc.
Ft. Lauderdale, Florida USA



CONTENT

Absorbers	4
BioSorb Soda Lime	
Anesthesia Mask	6
Standard Anesthesia Mask	
Scented Anesthesia Masks	
Airway Management	10
Endotracheal Tubes	
Endobronchial Tubes	
Endobronchial Blockers	
Laryngeal Masks	
Suction Catheter	
Closed Suction System	
Tracheostomy Tubes	
Endotracheal Tube Accessories	
Breathing Bags	42
Latex Breathing Bags	
Latex Free Breathing Bags	
Central Venous Catheter	44
Central Venous Catheter	
Anesthesia Spinal Needle	
Needles	48
Anesthesia Spinal Needle	
Epidural Needles	
Single use Epidural Needles	
Filters And Hme	54
Bacterial Viral Filters	
Heat Moisture Exchange Filter (HMEF)	
Procedural Trays	58
Anesthesia Procedural Trays	



Biotronix Healthcare®

BioSorb Soda Lime

Biotronix Healthcare's BioSorb Soda Lime is a carbon dioxide absorbent of choice for major hospitals worldwide. BioSorb is specifically manufactured for clinical use during anesthesia to absorb CO₂ within the breathing system. In these applications, breathable air is re-circulated for re-use and therefore any exhaled carbon dioxide needs to be removed in order to prevent a build-up to toxic levels.

Biosorb Soda Lime effectively removes carbon dioxide due to respiration carried out by the patient or the user.

Biosorb Benefits:

- Very low dust content
- High absorption capacity
- Color indicator white to violet
- As per appendix IX of MDD93/42/EEC
- Optimum hardness level 99% (As per USPXXII)

Biotronix Healthcare Biosorb Co₂ Absorber

REF #	DESCRIPTION	Quantity (per case)
AN10112-12	BioSorb CO ₂ Absorber 1.2kg Bottle	8Bot/Case
AN10112-45	BioSorb CO ₂ Absorber 4.5kg Bottle	2Bot/Case
AN10112-55	BioSorb CO ₂ Absorber 20kg Barrel	Each



Biotronix **Healthcare®**

Standard Anesthesia Mask

- Made from clear soft PVC
- Air filled cushion provide a comfortable fit
- Clear body provides easy patient observation
- Inflation valve for optimum seal
- Minimum pressure
- Latex free

Standar Anesthesia Mask

ITEM NO	SIZE	QTY/BOX
AN336001	Neonate	50
AN336002	Infant	50
AN336003	Child	50
AN336004	Large Child	50
AN336005	Small Adult	50
AN336006	Medium Adult	50
AN336007	Large Adult	50




Biotronix **Healthcare®**

Scented Anesthesia Masks

Biotronix Healthcare's disposable Scented Anesthesia Masks feature a clear, flexible tear drop design with a soft elastomeric cushion which provides an ergonomic fit and secure seal. During development, patient comfort and ease of use was our goal. Biotronix Healthcare's anesthesia masks come in neonate, toddler and child/small adult.

- Transparent body optimum clarity
- Profiled face cushion to facilitate secure seal on the patients face.
- Pre-inflated facial cushion manufactured from soft material.
- Top Mount Inflation valve for air pressure adjustment.
- Color Coded hook ring for use with head strap (Size 3)
- 22mm centrally located inlet.
- Scented masks come in bubble gum,



Biotronix Healthcare Scented Anesthesia Masks

REF #	DESCRIPTION	Quantity (per case)
AN337001	Scented Anesthesia Mask Neonate Size 1, Bubble Gum	50
AN337002	Scented Anesthesia Mask Neonate Size 1, Cherry	50
AN337003	Scented Anesthesia Mask Neonate Size 1, Strawberry	50
AN337004	Scented Anesthesia Mask Toddler Size 2, Bubble Gum	50
AN337005	Scented Anesthesia Mask Toddler Size 2, Cherry	50
AN337006	Scented Anesthesia Mask Toddler Size 2, Strawberry	50
AN337007	Scented Anesthesia Mask Child/Small Adult Size 3, Bubble Gum	50
AN337008	Scented Anesthesia Mask Child/Small Adult Size 3, Cherry	50
AN337010	Scented Anesthesia Mask Child/Small Adult Size 3, Strawberry	50



Biotronix Healthcare®

Endotracheal Tubes

Biotronix Healthcare intubating devices offer healthcare facilities with a complete line they can feel confident with. Our product range covers standard endotracheal tubes, preformed tubes, reinforced tubes, cuffed or uncuffed tubes so you can select the right device for each patient situation.

Standard Uncuffed ET Tubes

- Smooth rounded distal tip to facilitate insertion
- For oral and nasal intubation
- Depth lines to reference depth of intubation
- Provides the perfect balance of kink resistance flexibility and stiffness



Standard Uncuffed Endotracheal Tube

REF #	DESCRIPTION	QUANTITY (per box)	Quantity (per case)
AN211020U	Standard Uncuffed Endotracheal Tube, 2.0mm	10/Box	100 Boxes/Case
AN211025U	Standard Uncuffed Endotracheal Tube, 2.5mm	10/Box	100 Boxes/Case
AN211030U	Standard Uncuffed Endotracheal Tube, 3.0mm	10/Box	100 Boxes/Case
AN211035U	Standard Uncuffed Endotracheal Tube, 3.5mm	10/Box	100 Boxes/Case
AN211040U	Standard Uncuffed Endotracheal Tube, 4.0mm	10/Box	100 Boxes/Case
AN211045U	Standard Uncuffed Endotracheal Tube, 4.5mm	10/Box	100 Boxes/Case
AN211050U	Standard Uncuffed Endotracheal Tube, 5.0mm	10/Box	100 Boxes/Case
AN211055U	Standard Uncuffed Endotracheal Tube, 5.5mm	10/Box	100 Boxes/Case
AN211060U	Standard Uncuffed Endotracheal Tube, 6.0mm	10/Box	100 Boxes/Case
AN211065U	Standard Uncuffed Endotracheal Tube, 6.5mm	10/Box	100 Boxes/Case
AN211070U	Standard Uncuffed Endotracheal Tube, 7.0mm	10/Box	100 Boxes/Case
AN211075U	Standard Uncuffed Endotracheal Tube, 7.5mm	10/Box	100 Boxes/Case
AN211080U	Standard Uncuffed Endotracheal Tube, 8.0mm	10/Box	100 Boxes/Case
AN211085U	Standard Uncuffed Endotracheal Tube, 8.5mm	10/Box	100 Boxes/Case
AN211090U	Standard Uncuffed Endotracheal Tube, 9.0mm	10/Box	100 Boxes/Case
AN211095U	Standard Uncuffed Endotracheal Tube, 9.5mm	10/Box	100 Boxes/Case
AN2110100U	Standard Uncuffed Endotracheal Tube, 10.0mm	10/Box	100 Boxes/Case
AN2110105U	Standard Uncuffed Endotracheal Tube, 10.5mm	10/Box	100 Boxes/Case
AN2110110U	Standard Uncuffed Endotracheal Tube, 11.0mm	10/Box	100 Boxes/Case

High Volume Low Pressure Cuffed

- High volume, low pressure cuff
- Smooth rounded distal tip to facilitate insertion
- For oral and nasal intubation
- Depth lines to reference depth of intubation
- Provides the perfect balance of kink



High Volume Low Pressure Cuffed Endotracheal Tube

REF #	DESCRIPTION	QUANTITY (per box)	Quantity (per case)
AN211040H	Endotracheal Tube High Volume Low Pressure Cuffed, 4.0mm	10/Box	100 Boxes/Case
AN211045H	Endotracheal Tube High Volume Low Pressure Cuffed, 4.5mm	10/Box	100 Boxes/Case
AN211050H	Endotracheal Tube High Volume Low Pressure Cuffed, 5.0mm	10/Box	100 Boxes/Case
AN211055H	Endotracheal Tube High Volume Low Pressure Cuffed, 5.5mm	10/Box	100 Boxes/Case
AN211060H	Endotracheal Tube High Volume Low Pressure Cuffed, 6.0mm	10/Box	100 Boxes/Case
AN211065H	Endotracheal Tube High Volume Low Pressure Cuffed, 6.5mm	10/Box	100 Boxes/Case
AN211070H	Endotracheal Tube High Volume Low Pressure Cuffed, 7.0mm	10/Box	100 Boxes/Case
AN211075H	Endotracheal Tube High Volume Low Pressure Cuffed, 7.5mm	10/Box	100 Boxes/Case
AN211080H	Endotracheal Tube High Volume Low Pressure Cuffed, 8.0mm	10/Box	100 Boxes/Case
AN211085H	Endotracheal Tube High Volume Low Pressure Cuffed, 8.5mm	10/Box	100 Boxes/Case
AN211090H	Endotracheal Tube High Volume Low Pressure Cuffed, 9.0mm	10/Box	100 Boxes/Case
AN211095H	Endotracheal Tube High Volume Low Pressure Cuffed, 9.5mm	10/Box	100 Boxes/Case
AN2110100H	Endotracheal Tube High Volume Low Pressure Cuffed, 10.0mm	10/Box	100 Boxes/Case
AN2110105H	Endotracheal Tube High Volume Low Pressure Cuffed, 10.5mm	10/Box	100 Boxes/Case
AN2110110H	Endotracheal Tube High Volume Low Pressure Cuffed, 11.0mm	10/Box	100 Boxes/Case

Oral Preformed ET Tube Uncuffed

- Preformed design to reduce kinking
- Smooth rounded distal tip to facilitate insertion
- Depth lines to reference depth of intubation
- Provides the perfect balance of kink resistance flexibility and stiffness



Oral Preformed Endotracheal Tube Uncuffed

REF #	DESCRIPTION	QUANTITY (per box)	Quantity (per case)
AN212030U	Endotracheal Tube Oral Preformed Uncuffed, 3.0mm	10/Box	100 Boxes/Case
AN212035U	Endotracheal Tube Oral Preformed Uncuffed, 3.5mm	10/Box	100 Boxes/Case
AN212040U	Endotracheal Tube Oral Preformed Uncuffed, 4.0mm	10/Box	100 Boxes/Case
AN212045U	Endotracheal Tube Oral Preformed Uncuffed, 4.5mm	10/Box	100 Boxes/Case
AN212050U	Endotracheal Tube Oral Preformed Uncuffed, 5.0mm	10/Box	100 Boxes/Case
AN212055U	Endotracheal Tube Oral Preformed Uncuffed, 5.5mm	10/Box	100 Boxes/Case
AN212060U	Endotracheal Tube Oral Preformed Uncuffed, 6.0mm	10/Box	100 Boxes/Case
AN212065U	Endotracheal Tube Oral Preformed Uncuffed, 6.5mm	10/Box	100 Boxes/Case
AN212070U	Endotracheal Tube Oral Preformed Uncuffed, 7.0mm	10/Box	100 Boxes/Case
AN212075U	Endotracheal Tube Oral Preformed Uncuffed, 7.5mm	10/Box	100 Boxes/Case
AN212080U	Endotracheal Tube Oral Preformed Uncuffed, 8.0mm	10/Box	100 Boxes/Case
AN212085U	Endotracheal Tube Oral Preformed Uncuffed, 8.5mm	10/Box	100 Boxes/Case
AN212090U	Endotracheal Tube Oral Preformed Uncuffed, 9.0mm	10/Box	100 Boxes/Case
AN212095U	Endotracheal Tube Oral Preformed Uncuffed, 9.5mm	10/Box	100 Boxes/Case
AN2120100U	Endotracheal Tube Oral Preformed Uncuffed, 10.0mm	10/Box	100 Boxes/Case

Oral Preformed ET Tube Cuffed

- High volume, low pressure cuff
- Preformed design to reduce kinking
- High volume, low pressure cuff
- Smooth rounded distal tip to facilitate insertion
- Depth lines to reference depth of intubation



Oral Preformed High Volume Low Pressure Cuffed Endotracheal Tube

REF #	DESCRIPTION	QUANTITY (per box)	Quantity (per case)
AN212030H	Endotracheal Tube Oral Preformed High Volume Low Pressure Cuffed, 3.0mm	10/Box	100 Boxes/Case
AN212035H	Endotracheal Tube Oral Preformed High Volume Low Pressure Cuffed, 3.5mm	10/Box	100 Boxes/Case
AN212040H	Endotracheal Tube Oral Preformed High Volume Low Pressure Cuffed, 4.0mm	10/Box	100 Boxes/Case
AN212045H	Endotracheal Tube Oral Preformed High Volume Low Pressure Cuffed, 4.5mm	10/Box	100 Boxes/Case
AN212050H	Endotracheal Tube Oral Preformed High Volume Low Pressure Cuffed, 5.0mm	10/Box	100 Boxes/Case
AN212055H	Endotracheal Tube Oral Preformed High Volume Low Pressure Cuffed, 5.5mm	10/Box	100 Boxes/Case
AN212060H	Endotracheal Tube Oral Preformed High Volume Low Pressure Cuffed, 6.0mm	10/Box	100 Boxes/Case
AN212065H	Endotracheal Tube Oral Preformed High Volume Low Pressure Cuffed, 6.5mm	10/Box	100 Boxes/Case
AN212070H	Endotracheal Tube Oral Preformed High Volume Low Pressure Cuffed, 7.0mm	10/Box	100 Boxes/Case
AN212075H	Endotracheal Tube Oral Preformed High Volume Low Pressure Cuffed, 7.5mm	10/Box	100 Boxes/Case
AN212080H	Endotracheal Tube Oral Preformed High Volume Low Pressure Cuffed, 8.0mm	10/Box	100 Boxes/Case
AN212085H	Endotracheal Tube Oral Preformed High Volume Low Pressure Cuffed, 8.5mm	10/Box	100 Boxes/Case
AN212090H	Endotracheal Tube Oral Preformed High Volume Low Pressure Cuffed, 9.0mm	10/Box	100 Boxes/Case
AN212095H	Endotracheal Tube Oral Preformed High Volume Low Pressure Cuffed, 9.5mm	10/Box	100 Boxes/Case
AN2120100H	Endotracheal Tube Oral Preformed High Volume Low Pressure Cuffed, 10.0mm	10/Box	100 Boxes/Case

Nasal Preformed ET Tube Uncuffed

- Preformed design to reduce kinking
- Smooth rounded distal tip to facilitate insertion
- Depth lines to reference depth of intubation
- Provides the perfect balance of kink resistance flexibility and stiffness



Nasal Preformed Endotracheal Tube Uncuffed

REF #	DESCRIPTION	QUANTITY (per box)	Quantity (per case)
AN213030U	Endotracheal Tube Nasal Preformed Uncuffed, 3.0mm	10/Box	100 Boxes/Case
AN213035U	Endotracheal Tube Nasal Preformed Uncuffed, 3.5mm	10/Box	100 Boxes/Case
AN213040U	Endotracheal Tube Nasal Preformed Uncuffed, 4.0mm	10/Box	100 Boxes/Case
AN213045U	Endotracheal Tube Nasal Preformed Uncuffed, 4.5mm	10/Box	100 Boxes/Case
AN213050U	Endotracheal Tube Nasal Preformed Uncuffed, 5.0mm	10/Box	100 Boxes/Case
AN213055U	Endotracheal Tube Nasal Preformed Uncuffed, 5.5mm	10/Box	100 Boxes/Case
AN213060U	Endotracheal Tube Nasal Preformed Uncuffed, 6.0mm	10/Box	100 Boxes/Case
AN213065U	Endotracheal Tube Nasal Preformed Uncuffed, 6.5mm	10/Box	100 Boxes/Case
AN213070U	Endotracheal Tube Nasal Preformed Uncuffed, 7.0mm	10/Box	100 Boxes/Case
AN213075U	Endotracheal Tube Nasal Preformed Uncuffed, 7.5mm	10/Box	100 Boxes/Case
AN213080U	Endotracheal Tube Nasal Preformed Uncuffed, 8.0mm	10/Box	100 Boxes/Case
AN213085U	Endotracheal Tube Nasal Preformed Uncuffed, 8.5mm	10/Box	100 Boxes/Case
AN213090U	Endotracheal Tube Nasal Preformed Uncuffed, 9.0mm	10/Box	100 Boxes/Case
AN213095U	Endotracheal Tube Nasal Preformed Uncuffed, 9.5mm	10/Box	100 Boxes/Case
AN2130100U	Endotracheal Tube Nasal Preformed Uncuffed, 10.0mm	10/Box	100 Boxes/Case

Nasal Preformed ET Tube Cuffed

- High volume, low pressure cuff
- Preformed design to reduce kinking
- High volume, low pressure cuff
- Smooth rounded distal tip to facilitate insertion
- Depth lines to reference depth of intubation



Nasal Preformed High Volume Low Pressure Endotracheal Tube Cuffed

REF #	DESCRIPTION	QUANTITY (per box)	Quantity (per case)
AN213030H	Endotracheal Tube Nasal Preformed High Volume Low Pressure Cuffed, 3.0mm	10/Box	100 Boxes/Case
AN213035H	Endotracheal Tube Nasal Preformed High Volume Low Pressure Cuffed, 3.5mm	10/Box	100 Boxes/Case
AN213040H	Endotracheal Tube Nasal Preformed High Volume Low Pressure Cuffed, 4.0mm	10/Box	100 Boxes/Case
AN213045H	Endotracheal Tube Nasal Preformed High Volume Low Pressure Cuffed, 4.5mm	10/Box	100 Boxes/Case
AN213050H	Endotracheal Tube Nasal Preformed High Volume Low Pressure Cuffed, 5.0mm	10/Box	100 Boxes/Case
AN213055H	Endotracheal Tube Nasal Preformed High Volume Low Pressure Cuffed, 5.5mm	10/Box	100 Boxes/Case
AN213060H	Endotracheal Tube Nasal Preformed High Volume Low Pressure Cuffed, 6.0mm	10/Box	100 Boxes/Case
AN213065H	Endotracheal Tube Nasal Preformed High Volume Low Pressure Cuffed, 6.5mm	10/Box	100 Boxes/Case
AN213070H	Endotracheal Tube Nasal Preformed High Volume Low Pressure Cuffed, 7.0mm	10/Box	100 Boxes/Case
AN213075H	Endotracheal Tube Nasal Preformed High Volume Low Pressure Cuffed, 7.5mm	10/Box	100 Boxes/Case
AN213080H	Endotracheal Tube Nasal Preformed High Volume Low Pressure Cuffed, 8.0mm	10/Box	100 Boxes/Case
AN213085H	Endotracheal Tube Nasal Preformed High Volume Low Pressure Cuffed, 8.5mm	10/Box	100 Boxes/Case
AN213090H	Endotracheal Tube Nasal Preformed High Volume Low Pressure Cuffed, 9.0mm	10/Box	100 Boxes/Case
AN213095H	Endotracheal Tube Nasal Preformed High Volume Low Pressure Cuffed, 9.5mm	10/Box	100 Boxes/Case
AN2130100H	Endotracheal Tube Nasal Preformed High Volume Low Pressure Cuffed, 10.0mm	10/Box	100 Boxes/Case

Reinforced ET Tube Uncuffed

- Reinforced steel wire tube wall reducing kinking
- Smooth rounded distal tip to facilitate insertion
- For oral and nasal intubation
- Depth lines to reference depth of intubation



Reinforced Endotracheal Tube Uncuffed

REF #	DESCRIPTION	QUANTITY (per box)	Quantity (per case)
AN214030U	Reinforced Uncuffed Endotracheal Tube, 3.0mm	10/Box	100 Boxes/Case
AN214035U	Reinforced Uncuffed Endotracheal Tube, 3.5mm	10/Box	100 Boxes/Case
AN214040U	Reinforced Uncuffed Endotracheal Tube, 4.0mm	10/Box	100 Boxes/Case
AN214045U	Reinforced Uncuffed Endotracheal Tube, 4.5mm	10/Box	100 Boxes/Case
AN214050U	Reinforced Uncuffed Endotracheal Tube, 5.0mm	10/Box	100 Boxes/Case
AN214055U	Reinforced Uncuffed Endotracheal Tube, 5.5mm	10/Box	100 Boxes/Case
AN214060U	Reinforced Uncuffed Endotracheal Tube, 6.0mm	10/Box	100 Boxes/Case
AN214065U	Reinforced Uncuffed Endotracheal Tube, 6.5mm	10/Box	100 Boxes/Case
AN214070U	Reinforced Uncuffed Endotracheal Tube, 7.0mm	10/Box	100 Boxes/Case
AN214075U	Reinforced Uncuffed Endotracheal Tube, 7.5mm	10/Box	100 Boxes/Case
AN214080U	Reinforced Uncuffed Endotracheal Tube, 8.0mm	10/Box	100 Boxes/Case
AN214085U	Reinforced Uncuffed Endotracheal Tube, 8.5mm	10/Box	100 Boxes/Case
AN214090U	Reinforced Uncuffed Endotracheal Tube, 9.0mm	10/Box	100 Boxes/Case
AN214095U	Reinforced Uncuffed Endotracheal Tube, 9.5mm	10/Box	100 Boxes/Case
AN2140100U	Reinforced Uncuffed Endotracheal Tube, 10.0mm	10/Box	100 Boxes/Case

Reinforced ET Tube Cuffed

- Reinforced steel wire tube wall reducing kinking
- High volume, low pressure cuff
- Smooth rounded distal tip to facilitate insertion
- For oral and nasal intubation
- Depth lines to reference depth of intubation



Reinforced Endotracheal Tube Cuffed

REF #	DESCRIPTION	QUANTITY (per box)	Quantity (per case)
AN214030H	Reinforced Endotracheal Tube High Volume Low Pressure Cuffed, 3.0mm	10/Box	100 Boxes/Case
AN214035H	Reinforced Endotracheal Tube High Volume Low Pressure Cuffed, 3.5mm	10/Box	100 Boxes/Case
AN214040H	Reinforced Endotracheal Tube High Volume Low Pressure Cuffed, 4.0mm	10/Box	100 Boxes/Case
AN214045H	Reinforced Endotracheal Tube High Volume Low Pressure Cuffed, 4.5mm	10/Box	100 Boxes/Case
AN214050H	Reinforced Endotracheal Tube High Volume Low Pressure Cuffed, 5.0mm	10/Box	100 Boxes/Case
AN214055H	Reinforced Endotracheal Tube High Volume Low Pressure Cuffed, 5.5mm	10/Box	100 Boxes/Case
AN214060H	Reinforced Endotracheal Tube High Volume Low Pressure Cuffed, 6.0mm	10/Box	100 Boxes/Case
AN214065H	Reinforced Endotracheal Tube High Volume Low Pressure Cuffed, 6.5mm	10/Box	100 Boxes/Case
AN214070H	Reinforced Endotracheal Tube High Volume Low Pressure Cuffed, 7.0mm	10/Box	100 Boxes/Case
AN214075H	Reinforced Endotracheal Tube High Volume Low Pressure Cuffed, 7.5mm	10/Box	100 Boxes/Case
AN214080H	Reinforced Endotracheal Tube High Volume Low Pressure Cuffed, 8.0mm	10/Box	100 Boxes/Case
AN214085H	Reinforced Endotracheal Tube High Volume Low Pressure Cuffed, 8.5mm	10/Box	100 Boxes/Case
AN214090H	Reinforced Endotracheal Tube High Volume Low Pressure Cuffed, 9.0mm	10/Box	100 Boxes/Case
AN214095H	Reinforced Endotracheal Tube High Volume Low Pressure Cuffed, 9.5mm	10/Box	100 Boxes/Case
AN2140100H	Reinforced Endotracheal Tube High Volume Low Pressure Cuffed, 10.0mm	10/Box	100 Boxes/Case

Endotracheal Tube Cuffed with Evacuation Suction Lumen

- Evacuation lumen to reduce the risk of VAP
- High volume, low pressure cuff
- Smooth rounded distal tip to facilitate insertion
- For oral and nasal intubation
- Depth lines to reference depth of intubation



Endotracheal Tube With Suction Evacuation Lumen Cuffed


REF #	DESCRIPTION	QUANTITY (per box)	Quantity (per case)
AN215060H	Endotracheal Tube with Suction Evacuation Lumen, Cuffed, 6.0mm	10/Box	100 Boxes/Case
AN215065H	Endotracheal Tube with Suction Evacuation Lumen, Cuffed, 6.5mm	10/Box	100 Boxes/Case
AN215070H	Endotracheal Tube with Suction Evacuation Lumen, Cuffed, 7.0mm	10/Box	100 Boxes/Case
AN215075H	Endotracheal Tube with Suction Evacuation Lumen, Cuffed, 7.5mm	10/Box	100 Boxes/Case
AN215080H	Endotracheal Tube with Suction Evacuation Lumen, Cuffed, 8.0mm	10/Box	100 Boxes/Case
AN215085H	Endotracheal Tube with Suction Evacuation Lumen, Cuffed, 8.5mm	10/Box	100 Boxes/Case
AN215090H	Endotracheal Tube with Suction Evacuation Lumen, Cuffed, 9.0mm	10/Box	100 Boxes/Case

Biotronix Healthcare®

Endobronchial Tubes

Anesthesia clinicians and patients can count on Biotronix Healthcare's broad range of endobronchial tubes. We carry right sided and left sided options with multiple French sizes to meet all your clinical needs.

- Made from medical grade PVC.
- Equipped with switch connector and two low-friction suction catheters.
- Available in left-sided and right-sided.
- The cuffs, two lumens and pilot balloons are imprinted with Tracheal (clear) or Bronchial (blue) for easy identification.
- Smooth inner lumens can make the bronchoscopes easier to go through.
- Pre-loaded stylet can help tube maintain the shape for easier intubation.



Biotronix Healthcare Double Lumen Endobronchial Tubes

REF #	DESCRIPTION	Quantity (per Box)	Quantity (per case)
AN512001	Double Lumen Endobronchial Tubes Left-Sided 28FR	1/Box	10/Case
AN512002	Double Lumen Endobronchial Tubes Left-Sided 32FR	1/Box	10/Case
AN512003	Double Lumen Endobronchial Tubes Left-Sided 35FR	1/Box	10/Case
AN512004	Double Lumen Endobronchial Tubes Left-Sided 37FR	1/Box	10/Case
AN512005	Double Lumen Endobronchial Tubes Left-Sided 39FR	1/Box	10/Case
AN512006	Double Lumen Endobronchial Tubes Left-Sided 41FR	1/Box	10/Case
AN512007	Double Lumen Endobronchial Tubes Right-Sided 28FR	1/Box	10/Case
AN512008	Double Lumen Endobronchial Tubes Right-Sided 32FR	1/Box	10/Case
AN512009	Double Lumen Endobronchial Tubes Right-Sided 35FR	1/Box	10/Case
AN512010	Double Lumen Endobronchial Tubes Right-Sided 37FR	1/Box	10/Case
AN512011	Double Lumen Endobronchial Tubes Right-Sided 39FR	1/Box	10/Case
AN512012	Double Lumen Endobronchial Tubes Right-Sided 41FR	1/Box	10/Case

Double Lumen Endobronchial Tubes DEHP Free


REF #	DESCRIPTION		Quantity (per case)
AN513001	Double Lumen Endobronchial Tubes Left-Sided, DEHP Free, 28FR	1/Box	10/Case
AN513002	Double Lumen Endobronchial Tubes Left-Sided, DEHP Free, 32FR	1/Box	10/Case
AN513003	Double Lumen Endobronchial Tubes Left-Sided, DEHP Free, 35FR	1/Box	10/Case
AN513004	Double Lumen Endobronchial Tubes Left-Sided, DEHP Free, 37FR	1/Box	10/Case
AN513005	Double Lumen Endobronchial Tubes Left-Sided, DEHP Free, 39FR	1/Box	10/Case
AN513006	Double Lumen Endobronchial Tubes Left-Sided, DEHP Free, 41FR	1/Box	10/Case
AN513007	Double Lumen Endobronchial Tubes Right-Sided, DEHP Free, 28FR	1/Box	10/Case
AN513008	Double Lumen Endobronchial Tubes Right-Sided, DEHP Free, 32FR	1/Box	10/Case
AN513009	Double Lumen Endobronchial Tubes Right-Sided, DEHP Free, 35FR	1/Box	10/Case
AN513010	Double Lumen Endobronchial Tubes Right-Sided, DEHP Free, 37FR	1/Box	10/Case
AN513011	Double Lumen Endobronchial Tubes Right-Sided, DEHP Free, 39FR	1/Box	10/Case
AN513012	Double Lumen Endobronchial Tubes Right-Sided, DEHP Free, 41FR	1/Box	10/Case

Biotronix Healthcare®

Endobronchial Blockers

Endobronchial blocker tube is a device which can be inserted down a tracheal tube after tracheal intubation so as to block off the right or left main bronchus of the lungs in order to be able to achieve a controlled one-sided ventilation of the lungs in thoracic surgery. The Biotronix Healthcare Endobronchial Blocker has been specifically designed with anesthesiologist in mind to provide an innovative design for lung isolation.

- Comfortable and Soft
- Ultra-Thin medical grade PU
- Creates a superior seal
- Reduces Micro-Aspiration
- Twist Lock Design makes this easy to operate with secure locking
- Pediatric to Adult



Biotronix Healthcare Double Lumen Endobronchial Tubes

REF #	DESCRIPTION	Quantity (per Box)	Quantity (per case)
AN512001	Double Lumen Endobronchial Tubes Left-Sided 28FR	1/Box	10/Case
AN512002	Double Lumen Endobronchial Tubes Left-Sided 32FR	1/Box	10/Case
AN512003	Double Lumen Endobronchial Tubes Left-Sided 35FR	1/Box	10/Case
AN512004	Double Lumen Endobronchial Tubes Left-Sided 37FR	1/Box	10/Case
AN512005	Double Lumen Endobronchial Tubes Left-Sided 39FR	1/Box	10/Case
AN512006	Double Lumen Endobronchial Tubes Left-Sided 41FR	1/Box	10/Case
AN512007	Double Lumen Endobronchial Tubes Right-Sided 28FR	1/Box	10/Case
AN512008	Double Lumen Endobronchial Tubes Right-Sided 32FR	1/Box	10/Case
AN512009	Double Lumen Endobronchial Tubes Right-Sided 35FR	1/Box	10/Case
AN512010	Double Lumen Endobronchial Tubes Right-Sided 37FR	1/Box	10/Case
AN512011	Double Lumen Endobronchial Tubes Right-Sided 39FR	1/Box	10/Case
AN512012	Double Lumen Endobronchial Tubes Right-Sided 41FR	1/Box	10/Case

Double Lumen Endobronchial Tubes DEHP Free

REF #	DESCRIPTION		Quantity (per case)
AN513001	Double Lumen Endobronchial Tubes Left-Sided, DEHP Free, 28FR	1/Box	10/Case
AN513002	Double Lumen Endobronchial Tubes Left-Sided, DEHP Free, 32FR	1/Box	10/Case
AN513003	Double Lumen Endobronchial Tubes Left-Sided, DEHP Free, 35FR	1/Box	10/Case
AN513004	Double Lumen Endobronchial Tubes Left-Sided, DEHP Free, 37FR	1/Box	10/Case
AN513005	Double Lumen Endobronchial Tubes Left-Sided, DEHP Free, 39FR	1/Box	10/Case
AN513006	Double Lumen Endobronchial Tubes Left-Sided, DEHP Free, 41FR	1/Box	10/Case
AN513007	Double Lumen Endobronchial Tubes Right-Sided, DEHP Free, 28FR	1/Box	10/Case
AN513008	Double Lumen Endobronchial Tubes Right-Sided, DEHP Free, 32FR	1/Box	10/Case
AN513009	Double Lumen Endobronchial Tubes Right-Sided, DEHP Free, 35FR	1/Box	10/Case
AN513010	Double Lumen Endobronchial Tubes Right-Sided, DEHP Free, 37FR	1/Box	10/Case
AN513011	Double Lumen Endobronchial Tubes Right-Sided, DEHP Free, 39FR	1/Box	10/Case
AN513012	Double Lumen Endobronchial Tubes Right-Sided, DEHP Free, 41FR	1/Box	10/Case

Biotronix Healthcare®

Laryngeal Masks

Biotronix Healthcare carries a comprehensive range of laryngeal masks for routine hospital procedures and emergency anesthetic procedures. Our Laryngeal Masks are designed for patient comfort and safety. All our products are quality tested and meet all ISO international standards.

Double Lumen Silicon Laryngeal Mask

- Soft and smooth cuff for easy insertion and removal
- Designed with a gastric access channel
- Designed for patient comfort and safety
- Soft and smooth cuff for easy insertion



Double Lumen Silicon Laryngeal Mask

REF #	DESCRIPTION	QUANTITY (per box)	Quantity (per case)
AN244010	Double Lumen Silicon Laryngeal Mask #1.0, Gastric Tube 8FR	10/Box	4 Boxes/Case
AN244020	Double Lumen Silicon Laryngeal Mask #2.0, Gastric Tube 8FR	10/Box	4 Boxes/Case
AN244030	Double Lumen Silicon Laryngeal Mask #3.0, Gastric Tube 16FR	10/Box	4 Boxes/Case
AN244040	Double Lumen Silicon Laryngeal Mask #4.0, Gastric Tube 16FR	10/Box	4 Boxes/Case
AN244050	Double Lumen Silicon Laryngeal Mask #5.0, Gastric Tube 16FR	10/Box	4 Boxes/Case

Disposable Flexible Wire Reinforced Silicone Laryngeal Mask

- Suitable for adult children and infant use
- Stainless Steel Reinforced
- Wire reinforced tube reduces risk of kinking and airway tube occlusion
- Flexible tube allows for placement away from surgical field and improves surgical access while maintaining a tight seal
- Soft and smooth cuff for easy insertion and



Disposable Flexible Wire Reinforced Silicone Laryngeal Mask

REF #	DESCRIPTION	QUANTITY (per box)	Quantity (per case)
AN243010	Disposable Flexible Wire Reinforced Silicone Laryngeal Mask #1.0	10/Box	4 Boxes/Case
AN241015	Disposable Flexible Wire Reinforced Silicone Laryngeal Mask #1.5	10/Box	4 Boxes/Case
AN241020	Disposable Flexible Wire Reinforced Silicone Laryngeal Mask #2.0	10/Box	4 Boxes/Case
AN241025	Disposable Flexible Wire Reinforced Silicone Laryngeal Mask #2.5	10/Box	4 Boxes/Case
AN241030	Disposable Flexible Wire Reinforced Silicone Laryngeal Mask #3.0	10/Box	4 Boxes/Case
AN241040	Disposable Flexible Wire Reinforced Silicone Laryngeal Mask #4.0	10/Box	4 Boxes/Case
AN241050	Disposable Flexible Wire Reinforced Silicone Laryngeal Mask #5.0	10/Box	4 Boxes/Case

Disposable Silicone Laryngeal Mask

- Less traumatic than ET tube
- Convenient depth markers for correct positioning of the mask
- Soft and smooth cuff for easy insertion and removal
- Suitable for adult, children and infant use
- Latex free



Disposable Flexible Wire Reinforced Silicone Laryngeal Mask

REF #	DESCRIPTION	QUANTITY (per box)	Quantity (per case)
AN241010	Disposable Silicone Laryngeal Mask #1.0	10/Box	4 Boxes/Case
AN241015	Disposable Silicone Laryngeal Mask #1.5	10/Box	4 Boxes/Case
AN241020	Disposable Silicone Laryngeal Mask #2.0	10/Box	4 Boxes/Case
AN241025	Disposable Silicone Laryngeal Mask #2.5	10/Box	4 Boxes/Case
AN241030	Disposable Silicone Laryngeal Mask #3.0	10/Box	4 Boxes/Case
AN241040	Disposable Silicone Laryngeal Mask #4.0	10/Box	4 Boxes/Case
AN241050	Disposable Silicone Laryngeal Mask #5.0	10/Box	4 Boxes/Case

Reusable Silicone Laryngeal Mask

- Suitable for adult children and infant use
- Manufactured with high quality silicone
- Designed for patient comfort and safety
- Autoclavable
- Designed to be sterilized and used up to 40 times
- User card for patient information and usage count

Disposable Silicone Laryngeal Mask

REF #	DESCRIPTION	QUANTITY (per box)	Quantity (per case)
AN242010	Reusable Silicone Laryngeal Mask #1.0	10/Box	4 Boxes/Case
AN242015	Reusable Silicone Laryngeal Mask #1.5	10/Box	4 Boxes/Case
AN242020	Reusable Silicone Laryngeal Mask #2.0	10/Box	4 Boxes/Case
AN242025	Reusable Silicone Laryngeal Mask #2.5	10/Box	4 Boxes/Case
AN242030	Reusable Silicone Laryngeal Mask #3.0	10/Box	4 Boxes/Case
AN242040	Reusable Silicone Laryngeal Mask #4.0	10/Box	4 Boxes/Case
AN242050	Reusable Silicone Laryngeal Mask #5.0	10/Box	4 Boxes/Case
		10/Box	4 Boxes/Case



Biotronix Healthcare®

Suction Catheter

Suction catheters are intended for oropharyngeal, nasopharyngeal and tracheobronchial suctioning. Biotronix Healthcare's suction catheter products are designed to be soft and atraumatic for maximum patient comfort.

- Designed for tracheal/oropharyngeal/nasopharyngeal suction.
- Fingertip, thumb, and funnel control type connectors available for all your suctioning needs.
- Color-coded connectors at proximal end for instant identification of French catheter size.
- Smooth clear kink resistant PVC ensures efficient mucus suctioning while minimizing potential damage to the delicate tracheal mucous membranes.
- The large opening of the distal tip provides for easy and efficient removal of secretions.
- The 2 lateral eyes are smooth and close to the tip which maximizes suction efficiency.



Suction Catheter
Standard Connector



Suction Catheter
Fingertip Control



Suction Catheter
Thumb Control

Suction Catheter Standard Connector

REF #	DESCRIPTION	QUANTITY (per box)	Quantity (per case)
AN371005	Suction Catheter Standard Connector Size 5	100/Box	10 Boxes/Case
AN371006	Suction Catheter Standard Connector Size 6	100/Box	10 Boxes/Case
AN371008	Suction Catheter Standard Connector Size 8	100/Box	10 Boxes/Case
AN371010	Suction Catheter Standard Connector Size 10	100/Box	10 Boxes/Case
AN371012	Suction Catheter Standard Connector Size 12	100/Box	10 Boxes/Case
AN371014	Suction Catheter Standard Connector Size 14	100/Box	10 Boxes/Case
AN371016	Suction Catheter Standard Connector Size 16	100/Box	10 Boxes/Case
AN371018	Suction Catheter Standard Connector Size 18	100/Box	10 Boxes/Case


Suction Catheter Fingertip Control

REF #	DESCRIPTION	QUANTITY (per box)	Quantity (per case)
AN372005	Suction Catheter Fingertip Control Size 5	100/Box	10 Boxes/Case
AN372006	Suction Catheter Fingertip Control Size 6	100/Box	10 Boxes/Case
AN372008	Suction Catheter Fingertip Control Size 8	100/Box	10 Boxes/Case
AN372010	Suction Catheter Fingertip Control Size 10	100/Box	10 Boxes/Case
AN372012	Suction Catheter Fingertip Control Size 12	100/Box	10 Boxes/Case
AN372014	Suction Catheter Fingertip Control Size 14	100/Box	10 Boxes/Case
AN372016	Suction Catheter Fingertip Control Size 16	100/Box	10 Boxes/Case
AN372018	Suction Catheter Fingertip Control Size 18	100/Box	10 Boxes/Case



Suction Catheter Thumb Control

REF #	DESCRIPTION	QUANTITY (per box)	Quantity (per case)
AN373005	Suction Catheter Thumb Control Size 5	100/Box	10 Boxes/Case
AN373006	Suction Catheter Thumb Control Size 6	100/Box	10 Boxes/Case
AN373008	Suction Catheter Thumb Control Size 8	100/Box	10 Boxes/Case
AN373010	Suction Catheter Thumb Control Size 10	100/Box	10 Boxes/Case
AN373012	Suction Catheter Thumb Control Size 12	100/Box	10 Boxes/Case
AN373014	Suction Catheter Thumb Control Size 14	100/Box	10 Boxes/Case
AN373016	Suction Catheter Thumb Control Size 16	100/Box	10 Boxes/Case
AN373018	Suction Catheter Thumb Control Size 18	100/Box	10 Boxes/Case




Biotronix Healthcare®

Closed Suction System

Adult 24 Hour DOUBLE SWIVEL ELBOWS

- Gamma Sterilized
- Numbered Catheter Allowing the Clinician
- One-Way Lavage / Irrigation Port to Limit the Distance of Catheter Advance
- Dual Swivels Provides Flexibility Without Connectors
- Tracheostomy Sizes
- Tactile Sleeve for Better Feel
- With MDI PortClosed




Adult 24 Hour Double Swivel Elbows

PRODUCT CODE	OUTER DIAMETER	LENGTH	DESCRIPTION	PACKAGING
ANCS1010A	10 Fr 1 3.3mm	21.3" / 54cm	Double Swivel Elbow	25/CS
ANCS1010B	10 Fr 1 3.3mm	12" 130.5cm	Double Swivel Elbow Tracheostomy	25/CS
ANCS1012A	12 Fr 14mm	21.3" / 54cm	Double Swivel Elbow	25/CS
ANCS1012B	12Fr 14mm	12" 1 30.5cm	Double Swivel Elbow Tracheostomy	25/CS
ANCS1014A	14 Fr 1 4 .6mm	21.3" / 54cm	Double Swivel Elbow	25/CS
ANCS1014B	14 Fr 1 4 .6mm	12" 1 30.5cm	Double Swivel Elbow Tracheostomy	25/CS
ANCS1016A	16 Fr 1 5.3mm	21.3" / 54cm	Double Swivel Elbow	25/CS
ANCS1016B	16 Fr 1 5.3mm	12" 1 30.5cm	Double Swivel Elbow Tracheostomy	25/CS

Neonatal and Pediatric 24 Hour

- Gamma Sterilized
- Numbered Catheter Allowing the Clinician
- One-Way Lavage / Irrigation Port to Limit the Distance of Catheter Advance
- Tactile Sleeve for Better Feel



Neonatal and Pediatric 24 hour

PRODUCT CODE	OUTER DIAMETER	LENGTH	DESCRIPTION	PACKAGING
ANCS2005	5 Fr/,1.66mm	12" / 30.5cm	Y" Endotracheal,Tube Adapters (2.0mm, 2.5mm)	25/CS
ANCS2006	6 Fr/,2mm	12" / 30.5cm	Y" Endotracheal,Tube Adapters (2.5mm, 3.0mm, 3.5mm)	25/CS
ANCS2008	8 Fr/,2.6mm	12 "/,30.5cm	"Y" Endotracheal, Tube Adapters (3.5mm, 4.0mm, 4.5mm)	25/CS

Adult 72 hour Doble Swivel Elbows

- Gamma Sterilized
- Numbered Catheter Allowing the Clinician
- One-Way Lavage / Irrigation Port to Limit the Distance of Catheter Advance
- Tactile Sleeve for Better Feel

Adult 72 Hour Doble Swivel Elbows

PRODUCT CODE	OUTER DIAMETER	LENGTH	DESCRIPTION	PACKAGING
ANCS3010A	10 Fr/ 3.3mm	21.3" / 54cm	Double Swivel Elbow	25/CS
ANCS3010B	10 Fr/ 3.3mm	12" / 30.5cm	Double Swivel Elbow Tracheostomy	25/CS
ANCS3012A	12 Fr/ 4mm	21.3" / 54cm	Double Swivel Elbow	25/CS
ANCS3012B	12Fr / 4mm	12" / 30.5cm	Double Swivel Elbow Tracheostomy	25/CS
ANCS3014A	14 Fr/ 4.6mm	21.3" / 54cm	Double Swivel Elbow	25/CS
ANCS3014B	14 Fr/ 4.6mm	12" / 30.5cm	Double Swivel Elbow Tracheostomy	25/CS
ANCS3016A	16 Fr/ 5.3mm	21.3" / 54cm	Double Swivel Elbow	25/CS
ANCS3016B	16 Fr/ 5.3mm	12" / 30.5cm	Double Swivel Elbow Tracheostomy	25/CS



Biotronix
Healthcare®

Closed Suction System

Tracheostomy Tubes

Tracheostomy Tube is a hollow tube, with or without a cuff, that is electively inserted directly into the trachea through a surgical incision or with a wire-guided progressive dilatation technique in case of an emergency. Biotronix Healthcare Tracheostomy Tubes are designed in lightweight medical grade PVC for patient comfort.

Tracheostomy Tube-Uncuffed

- Designed in lightweight medical grade PVC for patient comfort
- Full length Radio-opaque line for correct placement detection
- Smooth rounded tip of the Obturator reduces trauma during insertion.
- Printed neck plate with size information for easy identification



Biotronix Healthcare Tracheostomy Tube Uncuffed

REF #	DESCRIPTION	Quantity (per Box)	Quantity (per case)
AN710030	Tracheostomy Tube Uncuffed, 3.0mm	1/Box	20/Case
AN710035	Tracheostomy Tube Uncuffed, 6.5mm	1/Box	20/Case
AN710040	Tracheostomy Tube Uncuffed, 4.0mm	1/Box	20/Case
AN710045	Tracheostomy Tube Uncuffed, 4.5mm	1/Box	20/Case
AN710050	Tracheostomy Tube Uncuffed 5.0mm	1/Box	20/Case
AN710055	Tracheostomy Tube Uncuffed, 5.5mm	1/Box	20/Case
AN710060	Tracheostomy Tube Uncuffed, 6.0mm	1/Box	20/Case
AN710065	Tracheostomy Tube Uncuffed, 6.5mm	1/Box	20/Case
AN710070	Tracheostomy Tube Uncuffed, 7.0mm	1/Box	20/Case
AN710075	Tracheostomy Tube Uncuffed, 7.5mm	1/Box	20/Case
AN710080	Tracheostomy Tube Uncuffed, 8.0mm	1/Box	20/Case
AN710085	Tracheostomy Tube Uncuffed, 8.5mm	1/Box	20/Case
AN710090	Tracheostomy Tube Uncuffed, 9.0mm	1/Box	20/Case
AN710095	Tracheostomy Tube Uncuffed, 9.5mm	1/Box	20/Case
AN710100	Tracheostomy Tube Uncuffed, 10.0mm	1/Box	20/Case

Tracheostomy Tube-High Volume Low Pressure Cuff

- Designed in lightweight medical grade PVC for patient comfort
- Full length Radio-opaque line for correct placement detection
- Smooth rounded tip of the Obturator reduces trauma during insertion.
- Printed neck plate with size information for easy identification



Biotronix Healthcare Tracheostomy Tube Cuffed

REF #	DESCRIPTION		Quantity (per case)
AN720030	Tracheostomy Tube High Volume Low Pressure Cuffed, 3.0mm	1/Box	20/Case
AN720035	Tracheostomy Tube High Volume Low Pressure Cuffed, 3.5mm	1/Box	20/Case
AN720040	Tracheostomy Tube High Volume Low Pressure Cuffed, 4.0mm	1/Box	20/Case
AN720045	Tracheostomy Tube High Volume Low Pressure Cuffed, 4.5mm	1/Box	20/Case
AN720050	Tracheostomy Tube High Volume Low Pressure Cuffed, 5.0mm	1/Box	20/Case
AN720055	Tracheostomy Tube High Volume Low Pressure Cuffed, 5.5mm	1/Box	20/Case
AN720060	Tracheostomy Tube High Volume Low Pressure Cuffed, 6.0mm	1/Box	20/Case
AN720065	Tracheostomy Tube High Volume Low Pressure Cuffed, 6.5mm	1/Box	20/Case
AN720070	Tracheostomy Tube High Volume Low Pressure Cuffed, 7.0mm	1/Box	20/Case
AN720075	Tracheostomy Tube High Volume Low Pressure Cuffed, 7.5mm	1/Box	20/Case
AN720080	Tracheostomy Tube High Volume Low Pressure Cuffed, 8.0mm	1/Box	20/Case
AN720085	Tracheostomy Tube High Volume Low Pressure Cuffed, 8.5mm	1/Box	20/Case
AN720090	Tracheostomy Tube High Volume Low Pressure Cuffed, 9.0mm	1/Box	20/Case
AN720095	Tracheostomy Tube High Volume Low Pressure Cuffed, 9.5mm	1/Box	20/Case
AN720100	Tracheostomy Tube High Volume Low Pressure Cuffed, 10.0mm	1/Box	20/Case

Tracheostomy Masks

- Clear, soft vinyl for patient comfort and visual assessment
- Helix Swivel 22mm Connector
- Secure Strap



Biotronix Healthcare Tracheostomy Masks

REF #	DESCRIPTION	Quantity (per case)
AN415001	Tracheostomy Mask, Adult,	50/Case
AN415006	Tracheostomy Mask, Pediatric,	50/Case



Endotracheal Tube Accessories

Biotronix Healthcare offers a variety of endotracheal tube accessories to meet all your patient needs. Comfort and ease of use is what we strive for.

Endotracheal Tube Holder

- Keep endotracheal tube in place without applying tape to the patients face.
- Accommodates multiple-size ET tubes.
- Bite block can protect patient teeth and tube from occlusion
- Allow suction of the oropharynx put into use.



Endotracheal Tube Holder

REF #	DESCRIPTION	Quantity (per box)	Quantity (per case)
ANS11002	Endotracheal Tube Holder, Medium	10/Box	100/Case
ANS11003	Endotracheal Tube Holder, Large	10/Box	100/Case

Tracheostomy Tube Holder

- Latex-free, foam-padded cotton for patient comfort.
- Velcro hooks attaching easily to neck flange.
- Wide band tapers at hooks to reduce excess secretion accumulation.
- Soft materials and wide band provide better skin care.



Tracheostomy Tube Holder

REF #	DESCRIPTION	Quantity (per box)	Quantity (per case)
ANS11004	Tracheostomy Tube Holder, Small	10/Box	100/Case
ANS11005	Tracheostomy Tube Holder, Medium	10/Box	100/Case
ANS11006	Tracheostomy Tube Holder, Large	10/Box	100/Case

ET Intubation Stylet

- Made of malleable aluminum covered with plastic sheath.
- Soft rounded tip for minimizing patient trauma during intubation.
- Flexible for assistance in intubation.
- Single use, sterile



ET Intubation Stylet

REF #	DESCRIPTION	Quantity (per box)	Quantity (per case)
ANS11008	ET Intubating Stylet 6FR	20/Box	100/Case
ANS11009	ET Intubating Stylet 10FR	20/Box	100/Case
ANS11010	ET Intubating Stylet 12FR	20/Box	100/Case
ANS11011	ET Intubating Stylet 14FR	20/Box	100/Case

Latex Breathing Bags

Biotronix Healthcare Breathing Bags include a complete range of latex non-conductive and latex-free bags in sizes 0.5 1.0, 2.0 and 3.0 liter with connector that meets International standards and provide tight fit. The multiple pleat design allows more positive control of manual ventilation with a secure fit and grip.

- Uniform inflation with a consistent feel
- All breathing bags are inspected to international standards
- Highest level of quality



Biotronix Healthcare Latex Breathing Bags

REF #	DESCRIPTION	Quantity (per case)
AN880543	0.5 Liter Latex Free Breathing Bags	50
AN880527	1 Liter Latex Free Breathing Bags	50
AN880526	2 Liter Latex Free Breathing Bags	50
AN880525	3 Liter Latex Free Breathing Bags	50

Biotronix Healthcare®

Latex Free Breathing Bags

Biotronix Healthcare Latex Free Breathing Bags include a complete range of latex non-conductive and latex-free bags in sizes 0.5 1.0, 2.0 and 3.0 liter with connector that meets International standards and provide tight fit. The multiple pleat design allows more positive control of manual ventilation with a secure fit and grip.

- Uniform inflation with a consistent feel
- All breathing bags are inspected to international standards
- Highest level of quality



Biotronix Healthcare Latex Free Breathing Bags

REF #	DESCRIPTION	Quantity (per case)
AN880521	0.5 Liter Latex Free Breathing Bags	50
AN880523	1 Liter Latex Free Breathing Bags	50
AN880522	2 Liter Latex Free Breathing Bags	50
AN880524	3 Liter Latex Free Breathing Bags	50

Biotronix Healthcare®

Central Venous Catheter

Biotronix Healthcare central venous catheters combine high quality manufacturing with procedure specific products, saving you time and creating cost efficiencies. Custom procedural trays are available; our account executive will work closely with you to create the perfect tray which meets your unique requirements.

Ask your Biotronix Healthcare account representative how best we can serve your facility's needs.

- Soft, flexible, radiopaque tips are less traumatic to the vessel
- Clear visibility of the catheter tip.
- Multiple sizes and lumen configurations available
- Attachable suture wing allows anchorage at puncture site






Movable Clamp avoids blood back flow



A round, soft hub reduces friction for patient comfort



Complete Kits: supplied with accessories for seldinger technique insertion. Catheter, vessel dilator, guidewire, introducer needle and insertion caps.



Soft polyurethane catheter and tip ensures trauma reduction to vessel walls.




Multiple introducer needles available syringe type or Y type.



Catheter and Guidewire with numbered markings to facilitate depth insertion.




Radiopaque catheter and tip facilitates vessel location.



Green Soft tip significantly reduces the risk of vascular trauma.

Basic Kits Contains: Catheter with clamps, J-Flex guidewire, Introducer Needle, Injection caps, 5cc Syringe, Vessel Dilator

ITEM NO.	TYPE	LUMEN SIZE	CATHETER LENGTH CM	QTY/B
ANCVCI-14-15	Single-lumen	14G	15	10
ANCVCI-14-20	Single-lumen	14G	20	10
ANCVCI-14-30	Single-lumen	14G	30	10
ANCVCI-16-15	Single-lumen	16G	15	10
ANCVCI-16-20	Single-lumen	16G	20	10
ANCVCI-16-30	Single-lumen	16G	30	10
ANCVCI-18-15	Single-lumen	18G	15	10
ANCVCI-18-20	Single-lumen	18G	20	10
ANCVCI-18-30	Single-lumen	18G	30	10
ANCVCI-20-13	Single-lumen	20G	13	10
ANCVCI-20-20	Single-lumen	20G	20	10
ANVCV2-40-05	Double-lumen	4F	5	10
ANVCV2-40-08	Double-lumen	4F	8	10
ANVCV2-40-13	Double-lumen	4F	13	10
ANVCV2-50-08	Double-lumen	SF	8	10
ANVCV2-50-13	Double-lumen	SF	13	10
ANVCV2-50-20	Double-lumen	SF	20	10
ANVCV2-70-15	Double-lumen	7F	15	10
ANVCV2- 70-20	Double-lumen	7F	20	10
ANVCV2- 70-30	Double-lumen	7F	30	10
ANVCV3-55-08	Triple-lumen	5.SF	8	10
ANVCV3-55-13	Triple-lumen	5.SF	13	10
ANVCV3- 70-15	Triple-lumen	7F	15	10
ANVCV3- 70-20	Triple-lumen	7F	20	10




Biotronix
Healthcare®

Central Venous Catheter

Anesthesia Spinal Needle

Single Use Quincke Spinal Needles

- Sharp bevel for low puncture force.
- Double bevel enables atraumatic puncture of the dura and minimizes the risk of post-spinal headache
- Transparent hub for quick detection of cerebrospinal fluid backflow and precise positioning of the needle tip.
- Mandrel grip is color-coded according to needle size.




Quincke Spinal Needles

Product Code	Description	Packaging
ANQ9000	Quincke Spinal Needle 18 Ga. x 3½ (90 mm)	"50/Bx 20Bx/CS"
ANQ9001	Quincke Spinal Needle 22 Ga. x 3½ (90 mm)	"50/Bx 20Bx/CS"
ANQ9002	Quincke Spinal Needle 22 Ga. x 1½ in. (40 mm)	"50/Bx 20Bx/CS"
ANQ9003	Quincke Spinal Needle 22 Ga. x 3½ in. (90 mm)	"50/Bx 20Bx/CS"
ANQ9004	Quincke Spinal Needle 22 Ga. x 4¾ in. (120 mm)	"50/Bx 20Bx/CS"
ANQ9005	Quincke Spinal Needle 22 Ga. x 5 in. (130 mm)	"50/Bx 20Bx/CS"
ANQ9006	Quincke Spinal Needle 22 Ga. x 7 in. (180 mm)	"50/Bx 20Bx/CS"
ANQ9007	Quincke Spinal Needle 25 Ga. x 3½ in. (90 mm)	"50/Bx 20Bx/CS"
ANQ9008	Quincke Spinal Needle 25 Ga. x 4¾ in. (120 mm)	"50/Bx 20Bx/CS"
ANQ9009	Quincke Spinal Needle 26 Ga. x 3½ in. (90 mm)	"50/Bx 20Bx/CS"
ANQ9010	Quincke Spinal Needle 27 Ga. x 3½ in. (90 mm)	"50/Bx 20Bx/CS"

Single Use Pencil Point Spinal Needles

- Innovative needle point design causes minimal trauma of the dura mater
- Low PDPH rates and excellent tactile feedback
- High flow rate and easy CSF visualization
- Precisely shaped and positioned side hole maintains adequate flow and reduces the chances of straddling the dura



Pencil Spinal Needles

Product Code	Description	Packaging
ANP9001	Pencil Point Spinal Needle 22 Ga. x 3½ (90 mm)	50/Bx 20Bx/ CS
ANP9002	Pencil Point Spinal Needle 24 Ga. x 4 in (100 mm)	50/Bx 20Bx/ CS
ANP9003	Pencil Point Spinal Needle 25 Ga. x 1 in. (25 mm)	50/Bx 20Bx/ CS
ANP9004	Pencil Point Spinal Needle 25 Ga. x 2 in. (50 mm)	50/Bx 20Bx/ CS
ANP9005	Pencil Point Spinal Needle 25 Ga. x 3½ in. (90 mm)	50/Bx 20Bx/ CS
ANP9006	Pencil Point Spinal Needle 25 Ga. x 5 in. (130 mm)	50/Bx 20Bx/ CS
ANP9007	Pencil Point Spinal Needle 25 Ga. x 6.1 in. (156 mm)	50/Bx 20Bx/ CS
ANP9008	Pencil Point Spinal Needle 27 Ga. x 3½ in. (90 mm)	50/Bx 20Bx/ CS
ANP9009	Pencil Point Spinal Needle 27 Ga. x 5 in. (127 mm)	50/Bx 20Bx/ CS

Biotronix *Healthcare®*

Epidural Needles

Single Use Epidural Needles

- Precision head has been rounded and bent to reduce the risk of spinal dura piercing
- Conveniently placed depth markings for easy visualization
- Attached wings provide a stable surface and easy to hold
- Ergonomic ridged hourglass shaped hub
- Polished inner bevel



Single Use Epidural Needles

Product Code	Description	Packaging
ANE9000	Tuohy Epidural Needle, 16 Ga. x 3¼ in. (80 mm), Clear Hub and Wings	50/Bx, 20Bx/CS
ANE9001	Tuohy Epidural Needle, 17 Ga. x 3½ in. (90 mm), Clear Hub and Wings	50/Bx, 20Bx/CS
ANE9002	Tuohy Epidural Needle, 17 Ga. x 4 in. (100 mm), Clear Hub and Wings	50/Bx, 20Bx/CS
ANE9003	Tuohy Epidural Needle, 17 Ga. x 4½ in. (115 mm), Clear Hub and Wings	50/Bx, 20Bx/CS
ANE9004	Tuohy Epidural Needle, 17 Ga. x 6 in. (150 mm), Clear Hub and Wings	50/Bx, 20Bx/CS
ANE9005	Tuohy Epidural Needle, 18 Ga. x 2 in. (50 mm), Clear Hub and Wings	50/Bx, 20Bx/CS
ANE9006	Tuohy Epidural Needle, 18 Ga. x 3½ in. (90 mm), Clear Hub and Wings	50/Bx, 20Bx/CS
ANE9007	Tuohy Epidural Needle, 18 Ga. x 4 in. (100 mm), Clear Hub and Wings	50/Bx, 20Bx/CS
ANE9008	Tuohy Epidural Needle, 18 Ga. x 4¾ in. (120 mm), Clear Hub and Wings	50/Bx, 20Bx/CS
ANE9009	Tuohy Epidural Needle, 18 Ga. x 6 in. (150 mm), Clear Hub and Wings	50/Bx, 20Bx/CS
ANE9010	Tuohy Epidural Needle, 20 Ga. x 2 in. (50 mm), Clear Hub and Wings	50/Bx, 20Bx/CS
ANE9011	Tuohy Epidural Needle, 20 Ga. x 3½ in. (90 mm), Clear Hub and Wings	50/Bx, 20Bx/CS
ANE9012	Tuohy Epidural Needle, 20 Ga. x 4 in. (100 mm), Clear Hub and Wings	50/Bx, 20Bx/CS
ANE9013	Tuohy Epidural Needle, 20 Ga. x 6 in. (150 mm), Clear Hub and Wings	50/Bx, 20Bx/CS
ANE9014	Tuohy Epidural Needle, 22 Ga. x 3.15 in. (80 mm), Clear Hub and Wings	50/Bx, 20Bx/CS
ANE9015	Tuohy Epidural Needle, 22 Ga. x 4 in. (100 mm), Clear Hub and Wings	50/Bx, 20Bx/CS

Single use Epidural Needles

The vast experience in the manufacturing of needles and state of the art production processes ensure that Biotronix Healthcare's Regional Anesthesia Needles and Trays are of the highest quality.

Regional Anesthesia helps to reduce pain during and after surgical procedures. Reduced pain in gynecological procedures, orthopedic surgical procedures, cardiovascular procedures, and most surgical procedures also reduces complications. Biotronix Healthcare offers a wide range

Combined Spinal & Epidural Block Needles:

- Blending the advantages of both Pencil Point Spinal Needles and Epidural needles
- Featuring the excellence of exact effect and significant muscle relaxing effect
- Delivering small amounts of anesthetic for local anesthesia
- Minimizes the effect post -dural puncture headaches

Combined Spinal & Epidural Block Needles:

Product Code	Description	Packaging
ANES9000	17 Ga. x 3½ in. Tuohy with 25 Ga. x 5 in. Pencil-Point Spinal Needle	50/Bx, 20Bx/CS
ANES9001	18 Ga. x 3½ in. Tuohy with 27 Ga. x 5½ in. Pencil-Point Spinal Needle	50/Bx, 20Bx/CS
ANES9002	18 Ga. x 3½ in. Tuohy with 27 Ga. x 5 in. Pencil-Point Spinal Needle	50/Bx, 20Bx/CS
ANES9003	18 Ga. x 3½ in. Tuohy with 27 Ga. x 5 in. Quincke Spinal Needle	50/Bx, 20Bx/CS
ANES9004	18 Ga. x 4½ in. Tuohy with 27 Ga. x 6 in. Pencil-Point Spinal Needle	50/Bx, 20Bx/CS



Bacterial Viral Filters

Biotronix Healthcare's Bacterial Viral Filter's are designed to be used with breathing and anesthesia systems, protecting the patient, hospital personnel and equipment from microbial contamination or particulate matter. Our Bacterial Viral Filters are lightweight, with a rounded ergonomic polypropylene housing ensuring patient safety and comfort.

- 99.99% cross-contamination efficiency
- Protects the subject against cross-infection
- Low dead space
- Lightweight
- Gas monitoring Luer Port

Parameters:

Description	Bacterial/Viral Filter (BV)
Humidity Output	N/A
Filtration Efficiency	"BFE 99.99%
VFE 99.999%"	
Resistance @ 30LPM	60PA
Dead Space	32ML
Tidal Volume Range	250-1500ML
Connections	22M/15F à 22F/15M
Gas Monitoring Luer Port with Retainer Strap	Yes
Weight	24 Grams

Biotronix Healthcare Bacterial Viral Filters

REF #	DESCRIPTION	Quantity (per case)
BA1001	Bacterial Viral Filter with Capnography Port 99.99%	50/Case



Biotronix Healthcare®

Heat Moisture Exchange Filter (HMEF)



Biotronix Healthcare Heat Moisture Exchange Filters are designed with low resistance to flow and high moisture output for ventilated and spontaneously breathing patients. Our high efficient filters help protect patients, healthcare workers and equipment from microbes in breathing circuits. We are always striving for patient safety and comfort.

HMEF V1

- Filtration Efficiency: BFE 99.999% VFE 99.999%
- Low breathing resistance
- Gas Sampling port
- Low dead space

Description	Heat & Moisture Exchange Filter (HMEF)
Humidity Output	31mg/H ₂ O@VT500ml
Filtration Efficiency	"BFE 99.999%
VFE 99.999%"	
Resistance @ 30LPM	180PA
Dead Space	45ML
Tidal Volume Range	150-1500ML
Connections	22M/15F à 22F/15M
Gas Sampling Luer Port with Retainer Strap	Yes
Weight	30 Grams

Biotronix Healthcare Heat Moisture Exchange Filter

REF #	DESCRIPTION	Quantity (per case)
BA1050	Heat Moisture Exchange Filter with Capnography Port 99.999%	50/Case

HMEF V2 Angled

- Filtration Efficiency: BFE 99.999% VFE 99.999%
- Low breathing resistance
- Gas Sampling port
- Low dead space



Description	Heat & Moisture Exchange Filter (HMEF)
Humidity Output	31mg/H ₂ O@VT500ml
Filtration Efficiency	BFE 99.999%
Resistance @ 30LPM	VFE 99.999%
Dead Space	180PA
Tidal Volume Range	45ML
Connections	150-1500ML
Gas Sampling Luer Port with Retainer Strap	22M/15F à 22F/15M
Weight	Yes

Biotronix Healthcare Angled Heat Moisture Exchange Filter

REF #	DESCRIPTION	Quantity (per case)
BA1060	Angled Heat Moisture Exchange Filter with Capnography Port 99.999%	50/Case

Anesthesia Procedural Trays

The use of Procedural Trays allows healthcare providers to work more efficiently and potentially deliver higher quality healthcare. They can save hospitals time before, during and after surgery by cutting back on preparation time and making waste removal easier. Procedural Packs drive up the quality of healthcare and assist healthcare providers to work more efficiently while cutting cost through better inventory management.

Biotronix Healthcare Anesthesia Procedural Trays come conveniently packed with all the components required for each procedure. Reducing procedural time frame and increasing patient procedural success increases patient safety and comfort

Tracheostomy Care Trays

- Towel/Underpad 18" x 13", paper/paper/poly
- Gloves Transparent, Med, Powderfree w/Wallet
- Trach Gauze Non-Woven, 4" x 4", 6 Ply
- Trach Brush 5.25"
- Pipe Cleaner 5.5"
- Cotton Tipped Applicator 5"



Biotronix Healthcare Tracheostomy Care Trays

REF #	DESCRIPTION	Quantity (per case)
RE40021	Tracheostomy Care Tray	20/Case

Suction Catheter Trays

- Chimney Suction Catheter w/ Graduated tubing
- PVS Gloves Transparent, M, Sleeved, 1 PR
- Paper Envelope
- Pop-up Solution Cap
- Paper Pouch



Biotronix Healthcare Suction Catheter Trays

REF #	DESCRIPTION	Quantity (per case)
AN4106	Chimney Suction Catheter 6FR, w/ Graduated Tubing	50/Case
AN4108	Chimney Suction Catheter 8FR, w/ Graduated Tubing	50/Case
AN4110	Chimney Suction Catheter 10FR, w/ Graduated Tubing	50/Case
AN4112	Chimney Suction Catheter 12FR, w/ Graduated Tubing	50/Case
AN4114	Chimney Suction Catheter 14FR, w/ Graduated Tubing	50/Case
AN4116	Chimney Suction Catheter 16FR, w/ Graduated Tubing	50/Case
AN4118	Chimney Suction Catheter 18FR, w/ Graduated Tubing	50/Case
AN4119	Thumb Control Suction Catheter 6FR, w/ Graduated Tubing	50/Case
AN4120	Thumb Control Suction Catheter 8FR, w/ Graduated Tubing	50/Case
AN4121	Thumb Control Suction Catheter 10FR, w/ Graduated Tubing	50/Case
AN4122	Thumb Control Suction Catheter 12FR, w/ Graduated Tubing	50/Case
AN4123	Thumb Control Suction Catheter 14FR, w/ Graduated Tubing	50/Case
AN4124	Thumb Control Suction Catheter 16FR, w/ Graduated Tubing	50/Case
AN4125	Thumb Control Suction Catheter 18FR, w/ Graduated Tubing	50/Case



Biotronix *Healthcare®*

Biotronix Healthcare Industries Inc.
Ft. Lauderdale, Florida USA