



## BIOP20 Coaxial Needle for BIOP800 and BIOP600



Specially designed for use with BIOP800  
and BIOP600 biopsy System.

### Features:

- Engineered compatibility with BIOP800 and BIOP600 biopsy system facilitates efficiency and accuracy.
- Lightweight and ultra thin wall construction.
- Create a clear path to work through when performing multiple biopsies in the same area.
- Sized just one gauge larger than corresponding core biopsy needles.
- Color-coded depth stopper conveniently match the compatible BIOP800 and BIOP600 biopsy system.

PRODUCT	GAUGE	LENGTH 10CM	LENGTH 16CM	LENGTH 20CM	PACKAGING
BIOP20 COAXILE BIOPSY NEEDLE	14	BIO81066	BIO81070	BIO81074	20/BX
	16	BIO81067	BIO81071	BIO81075	20/BX
	18	BIO81068	BIO81072	BIO81076	20/BX
	20	BIO81069	BIO81073	BIO81077	20/BX

Manufactured For:  
Biotronix Healthcare Industries Inc.  
Miami, Florida USA