

**Biotronix**  
Healthcare®

# SKIN STAPLER CATALOG



# Biotronix Healthcare®

## Skin Staplers

The Biotronix Healthcare Skin Staplers is a high-quality instrument with excellent performance for external skin closures. Our Skin Staplers come pre-loaded with 35 wide or regular medical grade stainless steel staples for use in a broad range of surgical procedures. The Biotronix Healthcare Skin Stapler will meet and exceed your healthcare facility's needs.

- Medical Grade Stainless Steel
- 35 Regular or Wide Staples
- See thru Window for Clear Visibility
- Audible Click
- Ergonomic Design for Easy Application



## Biotronix Healthcare Skin Staplers

REF #	DESCRIPTION	SIZE	STAPLE DIAMETER	STAPLE SIZE AFTER STAPLING	Quantity (per box)	Quantity (per case)
OR31035W	Skin Stapler Wide	Wide	0.6mm	7.2mm (width) × 4.9mm (height)	10/Box	8 Box/Case
OR310035R	Skin Stapler Regular	Regular	0.5mm	5.9mm (width) × 3.9mm (height)	10/Box	8 Box/Case







# ***Biotronix*** *Healthcare®*

Biotronix Healthcare Industries Inc.  
Ft. Lauderdale, Florida USA