


Biotronix
Healthcare®

TRADITIONAL
WOUND CARE
CATALOG



Biotronix Healthcare Industries Inc.
Ft. Lauderdale, Florida USA



CONTENT

Gauze Sponges	4
Sterile Surgical Gauze Sponges	
Non-Sterile Surgical Gauze Sponges	
Non-Woven All-Purpose Sponges	
Laparotomy Sponges	
X-Ray Detectable Gauze Sponges	
Cotton Filled Gauze Sponges	
Non-Woven Drain & IV Sponges Sterile	
Cotton Rolls	
Round Fluffy Sponges	
Round Stick Sponges	
Gauze/Bandage Rolls	18
Conforming Gauze Bandage Rolls	
Ultra-Gauze Bandage Rolls	
Elastic Bandages w/ Self Closures	
Elastic Bandages w/ Clips	
Cohesive Self Adhering Bandages	
Packing Strips	26
General Woundcare Trays	28
Wound Closures	31

Biotronix
Healthcare®

Biotronix Healthcare Industries Inc.
Ft. Lauderdale, Florida USA



Biotronix Healthcare® Gauze Sponges

Sterile Surgical Gauze Sponges

Biotronix Healthcare Gauze Sponges are 100% woven cotton and meet all the requirements your facility needs on weight and performance. Our Sponges can be used for several applications, are highly absorbent and offer exceptional clinical performance. Available in USP VII/ Basic and sterile/non-sterile. Latex free.

Features

- 100% cotton gauze
- Packed in 1's and 2's
- Packed in peel-down pouches for convenience or indexed trays of 10 for easy handling
- No loose edges exposed
- Highly absorbent
- Available in USP VII and Basic
- Sterile



Basic Sterile Surgical Gauze Sponges

REF #	DESCRIPTION	Quantity (per box)	Quantity (per case)
WOG401228	Sterile, Basic Gauze Sponge, 2" x 2", 8-ply , 2/pk	50 Pack/Box	30 Boxes /Case
WOG4013312	Sterile, Basic Gauze Sponge, 3" x 3", 12-ply , 2/pk	50 Pack/Box	24 Boxes /Case
WOG401448	Sterile, Basic Gauze Sponge, 4" x 4", 8-ply , 2/pk	50 Pack/Box	12 Boxes /Case
WOG4014412	Sterile, Basic Gauze Sponge, 4" x 4", 12-ply , 2/pk	25 Pack/Box	24 Boxes /Case



USP VII Sterile Surgical Gauze Sponges

REF #	DESCRIPTION	Quantity (per box)	Quantity (per case)
WOG401228P	Surgical Gauze Sponge USP VII Sterile 2's - 2"x 2" 8 Ply	50/Box	30 Boxes /Case
WOG4013312P	Surgical Gauze Sponge USP VII Sterile 2's - 3"x 3" 12 Ply	25/Box	48 Boxes /Case
WOG401448P	Surgical Gauze Sponge USP VII Sterile 2's - 4"x 4" 8 Ply	25/Box	24 Boxes /Case
WOG4014412P	Surgical Gauze Sponge USP VII Sterile 2's - 4"x 4" 12 Ply	25/Box	24 Boxes /Case
WOG4018412P	Surgical Gauze Sponge USP VII Sterile 2's - 8"x 4" 12 Ply	25/Box	20 Boxes /Case
WOG4104408T	Surgical Gauze Sponge USP VII Sterile, - 4"x 4" 8 Ply, (10 Indexed Tray)	10/Tray	128 Tray/Case
WOG4104412T	Surgical Gauze Sponge USP VII Sterile 10's/Tray - 4"x 4" 12 Ply	10/Tray	128 Trays /Case
WOG4104416T	Surgical Gauze Sponge USP VII Sterile, - 4"x 4" 16 Ply, (10 Indexed Tray)	10/Tray	128 Tray/Case
WOG4108412T	Surgical Gauze Sponge USP VII Sterile, - 8"x 4" 12 Ply, (10 Indexed Tray)	10/Tray	80 Tray/Case
WOG4002212	Surgical Gauze Sponges USP VII Sterile 1's - 2"x 2" 12 Ply	100/Box	24 Boxes /Case
WOG4003312	Surgical Gauze Sponges USP VII Sterile 1's - 3"x 3" 12 Ply	100/Box	24 Boxes /Case
WOG4004412	Surgical Gauze Sponges USP VII Sterile 1's - 4"x 4" 12 Ply	100/Box	12 Boxes /Case

Non-Sterile Surgical Gauze Sponges

Features

- 100% cotton gauze
- No loose edges exposed
- 200pcs per pack
- Available in USP VII and Basic
- Non-Sterile



USP VII Non-Sterile Surgical Gauze Sponges

REF #	DESCRIPTION	Quantity (per box)	Quantity (per case)
WOZ3102P	Surgical Gauze Sponge USP VII Non-Sterile, - 2"x 2" 8 Ply"	200/Pack	25 Pack /Case
WOZ3103P	Surgical Gauze Sponge USP VII Non-Sterile, - 2"x 2" 12 Ply"	200/Pack	40 Pack /Case
WOZ3104SP	Surgical Gauze Sponge USP VII Non-Sterile, - 3"x 3" 8 Ply"	200/Pack	20 Pack /Case
WOZ3104P	Surgical Gauze Sponge USP VII Non-Sterile, - 3"x 3" 12 Ply"	200/Pack	20 Pack /Case
WOZ3106P	Surgical Gauze Sponge USP VII Non-Sterile, - 4"x 4" 8 Ply"	200/Pack	20 Pack /Case
WOZ3102SP	Surgical Gauze Sponge USP VII Non-Sterile, - 4"x 4" 12 Ply"	200/Pack	10 Pack /Case
WOZ3103SP	Surgical Gauze Sponge USP VII Non-Sterile, - 4"x 4" 16 Ply"	200/Pack	10 Pack /Case
WOZ3108P	Surgical Gauze Sponge USP VII Non-Sterile, - 8"x 4" 12 Ply"	200/Pack	10 Pack /Case

Basic Non-Sterile Surgical Gauze Sponges

REF #	DESCRIPTION	Quantity (per box)	Quantity (per case)
WOZ3102	Non-Sterile, Basic Gauze Sponge, 2" x 2", 8-ply	200/Pack	25 Packs /Case
WOZ3103	Non-Sterile, Basic Gauze Sponge, 2" x 2", 12-ply	200/Pack	40 Packs /Case
WOZ3104S	Non-Sterile, Basic Gauze Sponge, 3" x 3", 8-ply	200/Pack	20 Packs /Case
WOZ3104	Non-Sterile, Basic Gauze Sponge, 3" x 3", 12-ply	200/Pack	20 Packs /Case
WOZ3106	Non-Sterile, Basic Gauze Sponge, 4" x 4", 8-ply 2	200/Pack	20 Packs /Case
WOZ3102S	Non-Sterile, Basic Gauze Sponge, 4" x 4", 12-ply	200/Pack	10 Packs /Case



Biotronix
Healthcare[®]

Gauze Sponges

Biotronix Healthcare® Gauze Sponges

Non-Woven All-Purpose Sponges

Biotronix Healthcare Non-Woven Sponges are an economical alternative to traditional woven gauze. These absorbent rayon/polyester blend sponges are ideal for general cleaning and prepping. They also perform well when applying ointments. Available in standard 30G or Premium 40G and sterile and non-sterile. Our sponges allow quick wicking and superior performance.

Features

- Non-Woven material
- Superior absorbency
- Less linting than standard gauze
- Sterile and Non-Sterile



Non-Woven All-Purpose Sponges Sterile

REF #	DESCRIPTION	Quantity (per box)	Quantity (per case)	
WOH3002024-2SP	Non-Woven Sponge Sterile 2's - 2"x 2" 4 Ply	40 Grams"	50/Box	30 Boxes /Case
WOH3003034-2SP	Non-Woven Sponge Sterile 2's - 3"x 3" 4 Ply	40 Grams"	50/Box	24 Boxes /Case
WOH3004044-2SP	Non-Woven Sponge Sterile 2's - 4"x 4" 4 Ply	40 Grams"	25/Box	24 Boxes /Case
WOH3002024-2S	Non-Woven Sponge Sterile 2's - 2"x 2" 4 Ply	30 Grams"	50/Box	30 Boxes /Case
WOH3003034-2S	Non-Woven Sponge Sterile 2's - 3"x 3" 4 Ply	30 Grams"	50/Box	24 Boxes /Case
WOH3004044-2S	Non-Woven Sponge Sterile 2's - 4"x 4" 4 Ply	30 Grams"	25/Box	24 Boxes /Case

Non-Woven All-Purpose Sponges Non-Sterile

REF #	DESCRIPTION	Quantity (per pack)	Quantity (per case)	
WOH3002024P	Non-Woven Sponge - 2"x 2" 4 Ply	40 Grams"	200/Pack	20 Pack /Case
WOH3003034P	Non-Woven Sponge - 3"x 3" 4 Ply	40 Grams"	200/Pack	20 Pack /Case
WOH3004044P	Non-Woven Sponge - 4"x 4" 4 Ply	40 Grams"	200/Pack	10 Pack /Case
WOH3002024	Non-Woven Sponge - 2"x 2" 4 Ply	30 Grams"	200/Pack	20 Pack /Case
WOH3003034	Non-Woven Sponge - 3"x 3" 4 Ply	30 Grams"	200/Pack	20 Pack /Case
WOH3004044	Non-Woven Sponge - 4"x 4" 4 Ply	30 Grams"	200/Pack	10 Pack /Case



Laparotomy Sponges

Biotronix Healthcare Laparotomy Sponges are manufactured of high quality cotton with a tight weave. Our sponges are ideal for use in surgical procedures. They are sterile, highly-absorbent, and cross-stitched for durability.

Features

- Pre-washed and sterile
- 100% high quality cotton with tight weave
- Cross stitched and super tight edges to reduce linting and loose threads
- X-ray detectable
- Banded in 5's for easy handling
- Latex Free

Laparotomy Sponges X-Ray Detectable

REF #	DESCRIPTION	Quantity (per box)	Quantity (per case)
WOF1212	Lap Sponge, Sterile, X-Ray Detectable, Pre-washed 12" x 12"	5/Pack	40 Pack /Case
WOF1818	Lap Sponge, Sterile, X-Ray Detectable, Pre-washed 18" x 18"	5/Pack	40 Pack /Case
WOF0836	Lap Sponge, Sterile, X-Ray Detectable, Pre-washed 8" x 36"	5/Pack	40 Pack /Case



Biotronix
Healthcare[®]
Gauze Sponges

Biotronix Healthcare®

Gauze Sponges

X-Ray Detectable Gauze Sponges

Biotronix Healthcare's X-Ray Detectable Gauze Sponges are economical gauze sponges ideal for use in surgical procedures. They are sterile, highly-absorbent, with cut edges folded in to keep loose threads from the wound site.

Features

- USP Type VII gauze, 100% woven-cotton
- Radiopaque thread for easy detection under x-ray
- Highly absorbent
- Fine mesh 100% cotton gauze, virtually lint-free
- Latex free
- Packed in convenient 10-count plastic trays
- Sterile



X-Ray Detectable Gauze Sponges USP VII, Sterile

REF #	DESCRIPTION	Quantity (per box)	Quantity (per case)	
WOX13010	X-Ray Detectable Surgical Gauze Sponge USP VII, Sterile - 4" x 4" x 16Ply	10/Tray"	10/Tray	128 Trays /Case
WOX13011	X-Ray Detectable Surgical Gauze Sponge USP VII, Sterile - 4" x 4" x 32Ply	10/Tray"	10/Box	30 Trays/Case
WOX13012	X-Ray Detectable Surgical Gauze Sponge USP VII, Sterile - 8" x 4" x 12Ply	10/Tray"	10/Box	48 Trays /Case
WOX13013	X-Ray Detectable Surgical Gauze Sponge USP VII, Sterile - 8" x 4" x 16Ply	10/Tray"	10/Box	48 Trays /Case
WOX13014	X-Ray Detectable Surgical Gauze Sponge USP VII, Sterile - 8" x 4" x 24Ply	10/Tray"	10/Box	50 Trays /Case



Cotton Filled Gauze Sponges

Biotronix Healthcare's Cotton-Filled Gauze Sponges are constructed with a gauze outer layer and a pure cotton inner layer. This inner layer provides extra absorbency and added patient comfort during use.

Features

- USP Type VII gauze, 100% woven-cotton
- Radiopaque thread for easy detection under x-ray
- Highly absorbent
- Fine mesh 100% cotton gauze, virtually lint-free
- Latex free
- Packed in convenient 10-count plastic trays
- Sterile



Cotton Filled Gauze Sponges

REF #	DESCRIPTION	Quantity (per box)	Quantity (per case)	
WOC11001	Surgical Gauze Sponge - 2"x 2" - Cotton Filled	Sterile"	200/Bag	25 Bags /Case
WOC11002	Surgical Gauze Sponge - 2"x 2" - Cotton Filled	Non-Sterile"	200/Bag	25 Bags /Case
WOC11003	Surgical Gauze Sponge 2" x 2" -Cottone Filled	Sterile 2/Pack"	50 Packs/Box	24 Boxes/Case

Non-Woven Drain & IV Sponges Sterile

Biotronix Healthcare's Non-woven drain and T sponges are manufactured from a rayon/polyester blend. They are precut to comfortably fit around IV sites, tracheostomies, catheters and tubes.

Features

- Made from a rayon/ polyester blend
- Less linting than regular gauze
- Designed to wick moisture away from the site
- Pre-cut T-slit provides a tight fit around catheters, tubes, I.V.'s and tracheotomies
- Also available in 100% woven cotton
- Latex free

Non-Woven Drain & IV Sponges Sterile

REF #	DESCRIPTION	Quantity (per box)	Quantity (per case)		
WOI11020	Sterile, I.V. Sponge, 2" x 2"	6-ply	2/Pack"	35 Packs /Box	20 Boxes /Case
WOI11021	Sterile, Drain Sponge, 4" x 4"	6-ply	2/Pack"	25 Packs /Box	12 Boxes /Case
WOI11022	Sterile, Drain Sponge, 100% Cotton, 4" x 4"	16-ply	1/Pack"	25 Packs /Box	12 Boxes /Case

Biotronix Healthcare®

Gauze Sponges

Cotton Rolls

Flexible cotton roll holds its shape with excellent performance for a number of dental procedures.

Features

- Uniform in size and shape
- Packed 50 per roll
- Soft and pliable, Medium, #2 rolls made from 100% cotton.
- Highly absorbent cotton fiber with superior performance.
- Maintains their shape during use.
- Latex and chlorine free

Cotton Rolls

REF #	DESCRIPTION	Quantity (per box)	Quantity (per case)
WOI11030	Cotton Rolls - #2 Medium, 50 Rolls/Wrap	2000/Box	10 Boxes /Case

Biotronix Healthcare®

Round Fluffy Sponges



Features:

- Soft and absorbent cotton sponge, provides extra protection to wounds
- For use with prepping/cleansing wound sites, pre and post-surgery
- Latex Free

Biotronix Healthcare Round Fluffy Sponges

REF #	DESCRIPTION	Quantity (per bag)	Quantity (per case)
WOQ24009	Non-Sterile Fluff Round Sponge , 3G	200/Bag	2000/Case
WOQ24010	Non-Sterile Super Fluff Round Sponge 29x23, 3.5G	200/Bag	2000/Case



Biotronix Healthcare®

Round Stick Sponges



Features:

- X-ray detectable
- Soft and absorbent cotton sponge, provides extra protection to wounds
- Used for deep sponging or prepping

Biotronix Healthcare Round Stick Sponges , X-Ray Detectable

REF #	DESCRIPTION	Quantity (per bag)	Quantity (per case)
WOX12011	Non Sterile, X Ray Detectable, 1" Large Bulk"	200/Bag	1000/Case
WOX12011XB	Non Sterile, X Ray Detectable, 1 1/2" X-Large Bulk"	200/Bag	1000/Case
WOX12011XS	Sterile, X Ray Detectable, 1 1/2" Large"	10/Bag	1000/Case

Biotronix Healthcare®

Gauze/Bandage Rolls

Conforming Gauze Bandage Rolls

Biotronix Healthcare's gauze products address the needs of both patient and physicians. Biotronix is involved in the manufacturing and sales of over 20 types with over 200 specifications of gauze sponges, gauze rolls, dressings, OR Towels and bandages. Biotronix Healthcare is sure to be able to supply all your gauze product needs. Our commitment to high quality and strict compliance with international standards assures you of a quality product at an affordable price.

Features

Gently secures to any body contour without restriction of movement and accommodate an increase or decrease in swelling.

Uniquely constructed weave and finished edges eliminating unraveling and lint.

The combination of rayon and polyester help the bandage to provide stretch memory when applied.





Conforming Stretch Gauze Bandage Rolls

REF #	DESCRIPTION	Quantity (per box)	Quantity (per case)		
WO3102	Conforming Stretch Gauze Bandage Rolls 2" x 4.1yds	Non-Sterile	No Pkg"	12/Box	8 Box /Case
WO3103	Conforming Stretch Gauze Bandage Rolls 3" x 4.1yds	Non-Sterile	No Pkg"	12/Box	8 Box /Case
WO3104	Conforming Stretch Gauze Bandage Rolls 4" x 4.1yds	Non-Sterile	No Pkg"	12/Box	8 Box /Case
WO3106	Conforming Stretch Gauze Bandage Rolls 6" x 4.1yds	Non-Sterile	No Pkg"	6/Box	8 Box /Case
WO3102P	Conforming Stretch Gauze Bandage Rolls 2" x 4.1yds	Non-Sterile	Individually Packed*	12/Box	8 Box /Case
WO3103P	Conforming Stretch Gauze Bandage Rolls 3" x 4.1yds	Non-Sterile	Individually Packed*	12/Box	8 Box /Case
WO3104P	Conforming Stretch Gauze Bandage Rolls 4" x 4.1yds	Non-Sterile	Individually Packed*	12/Box	8 Box /Case
WO3106P	Conforming Stretch Gauze Bandage Rolls 6" x 4.1yds	Non-Sterile	Individually Packed*	6/Box	8 Box /Case
WO3102S	Conforming Stretch Gauze Bandage Rolls 2" x 4.1yds	Sterile	Peel Open Blister Pouch"	12/Box	8 Box /Case
WO3103S	Conforming Stretch Gauze Bandage Rolls 3" x 4.1yds	Sterile	Peel Open Blister Pouch"	12/Box	8 Box /Case
WO3104S	Conforming Stretch Gauze Bandage Rolls 4" x 4.1yds	Sterile	Peel Open Blister Pouch"	12/Box	8 Box /Case
WO3106S	Conforming Stretch Gauze Bandage Rolls 6" x 4.1yds	Sterile	Peel Open Blister Pouch"	6/Box	8 Box /Case

Ultra-Gauze Bandage Rolls



Features:

Ideal for primary and secondary wound care dressings.

When used as a primary dressing, the open weave design provides fast wicking, absorbency and aeration.

Has a tightly finished edge for reduced lint; and a unique crimp design for increased cushioning and patient comfort.

The 6-ply bandage stretches and conforms to a variety of patient needs.

Ultra Gauze Bandage Rolls

REF #	DESCRIPTION	Quantity (per box)	Quantity (per case)	
WO3161	Ultra Gauze Bandage Rolls, 4.5" x 4.1yds	6 Ply	Sterile"	100 /Case
WO3162	Ultra Gauze Bandage Rolls, 4.5" x 4.1yds	6 Ply	Non-Sterile"	100 /Case

Biotronix Healthcare®

Gauze/Bandage Rolls

Elastic Bandages w/ Self Closures

Features:

- Easy-to-use, fast and apply maximum compression
- Comfortable and strong enough for laundering and repeated usage
- Individually wrapped
- Easy use self-closure
- 4.5 yds stretched
- Latex free




Elastic Banadages w/ Self Closure

REF #	DESCRIPTION	Quantity (per box)	Quantity (per case)
WO5102	Elastic Banadages w/ Self Closure, 2" X 5yds"	5/Bag	10 Bag /Case
WO5103	Elastic Banadages w/ Self Closure, 3" X 5yds"	5/Bag	10 Bag /Case
WO5104	Elastic Banadages w/ Self Closure, 4" X 5yds"	5/Bag	10 Bag /Case
WO5106	Elastic Banadages w/ Self Closure, 6" X 5yds"	5/Bag	10 Bag/Case

Biotronix Healthcare®

Gauze/Bandage Rolls

Elastic Bandages w/ Clips



Features:

- Easy-to-use, fast and apply maximum compression
- Comfortable and strong enough for laundering and repeated usage
- Easy secure clips
- Individually wrapped
- Latex free

Elastic Bandage w/ Clips

REF #	DESCRIPTION	Quantity (per box)	Quantity (per case)
WO5202	Elastic Bandage w/ Clips, 2" X 5yds"	5/Bag	10 Bag /Case
WO5203	Elastic Bandage w/ Clips, 3" X 5yds"	5/Bag	10 Bag /Case
WO5204	Elastic Bandage w/ Clips, 4" X 5yds"	5/Bag	10 Bag /Case
WO5206	Elastic Bandage w/ Clips, 6" X 5yds"	5/Bag	10 Bag/Case



Cohesive Self Adhering Bandages

Features:

- Self-Adhering bandage
- Soft foam texture.
- Acts as a secondary bandage without the need for clips or tape
- Available in various colors
- Hand-tears to desired length for easy-application



Cohesive Self Adhering Bandages

REF #	DESCRIPTION	Quantity (per box)	Quantity (per case)
W08101	Cohesive Self Adhering Bandages 1" x 5yds. beige	30	30/Case
W08102	Cohesive Self Adhering Bandages 2" x 5yds. beige	30	30/Case
W08103	Cohesive Self Adhering Bandages 3" x 5yds. beige	24	24/Case
W08104	Cohesive Self Adhering Bandages 4" x 5yds. beige	18	18/Case

REF #	DESCRIPTION	Quantity (per box)	Quantity (per case)
WO8105	Cohesive Self Adhering Bandages 5" x 5yds. Beige*		12/Case
WO8106	Cohesive Self Adhering Bandages 1" x 5yds. Dark Blue*		30/Case
WO8107	Cohesive Self Adhering Bandages 2" x 5yds. Dark Blue*		38/Case
WO8108	Cohesive Self Adhering Bandages 3" x 5yds. Dark Blue*		24/Case
WO8109	Cohesive Self Adhering Bandages 4" x 5yds. Dark Blue*		18/Case
WO8110	Cohesive Self Adhering Bandages 5" x 5yds. Dark Blue*		12/Case
WO8111	Cohesive Self Adhering Bandages 1" x 5yds. Light Blue*		30/Case
WO8112	Cohesive Self Adhering Bandages 2" x 5yds. Light Blue*		38/Case
WO8113	Cohesive Self Adhering Bandages 3" x 5yds. Light Blue*		24/Case
WO8114	Cohesive Self Adhering Bandages 4" x 5yds. Light Blue*		18/Case
WO8115	Cohesive Self Adhering Bandages 5" x 5yds. Light Blue*		12/Case
WO8116	Cohesive Self Adhering Bandages 1" x 5yds. Red*		30/Case
WO8117	Cohesive Self Adhering Bandages 2" x 5yds. Red*		38/Case
WO8118	Cohesive Self Adhering Bandages 3" x 5yds. Red*		24/Case
WO8119	Cohesive Self Adhering Bandages 4" x 5yds. Red*		18/Case
WO8120	Cohesive Self Adhering Bandages 5" x 5yds. Red*		12/Case
WO8121	Cohesive Self Adhering Bandages 1" x 5yds. Purple*		30/Case
WO8122	Cohesive Self Adhering Bandages 2" x 5yds. Purple*		38/Case
WO8123	Cohesive Self Adhering Bandages 3" x 5yds. Purple*		24/Case
WO8124	Cohesive Self Adhering Bandages 4" x 5yds. Purple*		18/Case
WO8125	Cohesive Self Adhering Bandages 5" x 5yds. Purple*		12/Case
WO8126	Cohesive Self Adhering Bandages 1" x 5yds. Pink*		30/Case
WO8127	Cohesive Self Adhering Bandages 2" x 5yds. Pink*		38/Case



REF #	DESCRIPTION	Quantity (per box)	Quantity (per case)
WO8126	Cohesive Self Adhering Bandages 3" x 6yds Pink [†]		24/Case
WO8129	Cohesive Self Adhering Bandages 4" x 6yds Pink [†]		18/Case
WO8130	Cohesive Self Adhering Bandages 5" x 6yds Pink [†]		12/Case
WO8131	Cohesive Self Adhering Bandages 1" x 6yds Green [†]		30/Case
WO8132	Cohesive Self Adhering Bandages 2" x 6yds Green [†]		36/Case
WO8133	Cohesive Self Adhering Bandages 3" x 6yds Green [†]		24/Case
WO8134	Cohesive Self Adhering Bandages 4" x 6yds Green [†]		18/Case
WO8135	Cohesive Self Adhering Bandages 5" x 6yds Green [†]		12/Case

Biotronix Healthcare® Packing Strips

Packing Strips

Packing strips from Biotronix Healthcare are a sterile single-use wound dressing comprising in a single cotton gauze strip, packed in amber colored jars either plain or impregnated with iodoform solution. Their primary use is for drainage of open and/or infected wounds. Our tightly woven gauze strips will not lint or fray. The Biotronix Healthcare Packing Strips are used in operating rooms, surgical centers, oral surgery, nursing homes and emergency rooms.

Iodoform Packing Strip

- All natural 100% Cotton
- Fine mesh gauze ideal for wet and dry packing
- Impregnated with Iodoform solution
- All are packaged one strip per bottle, sterile
Nasal Packinga

Plain Packing Strips

- All natural 100% Cotton
- Fine mesh gauze ideal for wet and dry packing
- Packed plain
- All are packaged one strip per bottle, sterile
Nasal Packing
- Dental Packing





Packing Strips

REF #	DESCRIPTION	Quantity (per case)	
WO60-2505	Plain Packing Strip, 1/4" x 5yds	1/Bottle"	12/Case
WO60-5005	Plain Packing Strip, 1/2" x 5yds	1/Bottle"	12/Case
WO60-1005	Plain Packing Strip, 1" x 5yds	1/Bottle"	12/Case
WO60-2005	Plain Packing Strip, 2" x 5yds	1/Bottle"	12/Case
WO61-2505	Iodoform Packing Strip, 1/4" x 5yds	1/Bottle"	12/Case
WO61-5005	Iodoform Packing Strip, 1/2" x 5yds	1/Bottle"	12/Case
WO61-1005	Iodoform Packing Strip, 1" x 5yds	1/Bottle"	12/Case
WO61-2005	Iodoform Packing Strip, 2" x 5yds	1/Bottle"	12/Case

Biotronix Healthcare®

General Woundcare Trays

Woundcare Trays

The use of Procedural Trays allows healthcare providers to work more efficiently and potentially deliver higher quality healthcare. They can save hospitals time before, during and after surgery by cutting back on preparation time and making waste removal easier. Procedural Packs drive up the quality of healthcare and assist healthcare providers to work more efficiently while cutting cost through better inventory management.

Biotronix Healthcare General Hospital Trays come conveniently packed with all the components required for each procedure. Reducing procedural time frame and increasing patient procedural success increases patient safety and comfort.

Features & Benefits:


Minor Laceration Tray

- Compartmentalized Tray w/ Lid Towle
- Towel / Underpad 26" x 16", Blue
- Fenestrated Drape 26" x 16", White
- Absorbent Towel 18" x 17", White
- Gauze Sponge 4" x 4", 8 Ply
- Needle Holder Floor Grade, 5"
- Tissue Forceps with Teeth Wire Grade, 5"




Biotronix Healthcare Minor Laceration Tray

REF #	DESCRIPTION	Quantity (per case)
OR4004	Minor Laceration Tray	20/Case



Biotronix
Healthcare®

General Woundcare Trays




Wound Closure Devices

The Biotronix Healthcare AcuLok I and AcuLok II skin closure line is used to close surface wounds and surgical incisions in the final stage. Our devices can be used in conjunction with sutures or staple closure of the underlying tissue to close the skin layer. Our devices furnish a quick and superior result promoting healthy wound healing and avoiding infections and potential complications.

AcuLok II Wound Closure Device

- Quick application
- Eliminates Track marks for cosmetic appearance
- May eliminates follow up visit for removal
- Enhanced patient comfort



AcuLok II Wound Closure Clip Device

REF #	DESCRIPTION	Quantity (per envelope)	Quantity (per box)	Quantity (per case)
WO51000	AcuLok II Wound Closure Clip Device Regular	5 Clip/Envelope	12 Envelope/Box	8 Boxes/Case
WO51001	AcuLok II Wound Closure Clip Device Large	3 Clip/Envelope	12 Envelope/Box	8 Boxes/Case

Biotronix Healthcare®

Wound Closures

AcuLok I Wound Closure Strips

- Quick application
- Comfortable to wear
- Less Scarring
- Less chance of infection



AcuLok I Wound Closure Strips

REF #	DESCRIPTION	Quantity (per envelope)	Quantity (per box)	Quantity (per case)
WO52000	AcuLok I Wound Closure Strips 1/8" x 1-1/2"	8/Envelope	50/Box	4 Box/Case
WO52001	AcuLok I Wound Closure Strips 1/8" x 3"	5/Envelope	50/Box	4 Box/Case
WO52002	AcuLok I Wound Closure Strips 1/4" x 3"	3/Envelope	50/Box	4 Box/Case
WO52003	AcuLok I Wound Closure Strips 1/4" x 4"	10/Envelope	50/Box	4 Box/Case
WO52004	AcuLok I Wound Closure Strips 1/2" x 4"	6/Envelope	50/Box	4 Box/Case



Biotronix ***Healthcare®***

Biotronix Healthcare Industries Inc.
Ft. Lauderdale, Florida USA